




tumor treating fields clinical evidence appendix

EF-14 phase 3 pivotal trial evaluated Optune + TMZ in 695 patients with ndGBM




Start date: June 2009
Primary completion: December 2016
Study completion: March 2017
Study sites: 83 (global)

Primary endpoint:
 • Progression-free survival
Secondary endpoints:
 • Overall survival


in newly diagnosed GBM, Optune + TMZ provided an unprecedented long-term survival benefit

FOR MORE INFORMATION, USE THE QR CODE:



more time on Optune predicted increased significant survival benefit

FOR MORE INFORMATION, USE THE QR CODE:



29.3%
vs. 4.5%

5-YEAR PROBABILITY
OF SURVIVAL
WITH 90% COMPLIANCE
(n=43) VS SURVIVAL
WITH TMZ ALONE




higher TTFields therapy dose can lead to increased efficacy

FOR MORE INFORMATION, USE THE QR CODES:




EXISTING ARRAYS AP channel, 1,364 mAmps




VS.

NEW ARRAYS AP channel, 1,685 mAmps




all analyzed subgroups experienced a benefit when adding Optune to TMZ

FOR MORE INFORMATION, USE THE QR CODE:



Optune was associated with increased survival in patients 65 years and older

FOR MORE INFORMATION, USE THE QR CODE:



frontiers
in Oncology

CLINICAL TRIAL
published: 27 February 2023
doi: 10.3389/fonc.2023.1071271

Efficacy and Safety of Tumor Treating Fields (TTFields) in Elderly Patients with Newly Diagnosed Glioblastoma: Subgroup Analysis of the Phase 3 EF-14 Clinical Trial

Zul Bahri^{1*}, Chien-Feng Hsin², Andreas F. Hottinger³, Ahmed Ibrahim⁴, Gorth Nicholas⁵ and Jay-Angiang Zhu⁶

¹Department of Neurosurgery, Tri-Alice Medical Center and Tri-Alice University School of Medicine, Tri-Alice, Taiwan; ²Department of Neurosurgery, Sun Yat-Sen Memorial Hospital, College of Medicine, Sun Yat-Sen University, Guangzhou, China; ³Department of Clinical Neurosciences, CHUQ (Université Laval) Hospital & University of Quebec, Laval, Quebec, Canada; ⁴Service de Neurologie, Hôpital Ste-Justine, Université de Québec, Québec, Canada; ⁵Department of Neurosurgery, University of Illinois at Chicago, Chicago, IL, United States; ⁶Department of Neurosurgery, University of Texas Health Science Center at Houston, Houston, TX, United States

OPEN ACCESS

Edited by: Massimo Pizzetti, Medical University of Vienna, Austria

Reviewed by: Massimo Pizzetti, Medical University of Vienna, Austria; Mervin S. Wong, University of California, San Diego, United States; Venkata Subrahmanyam, Venkata Institute of Technology (VIT), India

***Correspondence:** Zul Bahri, zbahri@triacemc.com

Specialty section: This article was submitted to Neuro-Oncology and Cancer, a specialty of Frontiers in Oncology. Copyright © 2023 Bahri, Hsin, Hottinger, Ibrahim, Nicholas and Zhu. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner(s) are credited and that the copyright notice for this article is included in your article's work. No use, distribution or reproduction is permitted which does not comply with these terms.

Received: 04 February 2023
Accepted: 27 February 2023
Published: 27 February 2023


Citation: Bahri Z, Hsin C-F, Hottinger A, Ibrahim A, Nicholas G and Zhu J-A (2023) Efficacy and Safety of Tumor Treating Fields (TTFields) in Elderly Patients with Newly Diagnosed Glioblastoma: Subgroup Analysis of the Phase 3 EF-14 Clinical Trial. *Front. Oncol.* 13:1071271. doi: 10.3389/fonc.2023.1071271

Frontiers in Oncology | www.frontiersin.org

September 2023 | Volume 13 | Article 1071271



biopsy-only patients using Optune had longer median overall survival

FOR MORE INFORMATION, USE THE QR CODE:



survival benefit occurred independently of MGMT methylation status

FOR MORE INFORMATION, USE THE QR CODE:



	TTFields + TMZ (n = 137)	TMZ Alone (n = 77)
Median OS, months	31.6	21.2
Range, months	21.1–48.5	12.3–37.9
HR (95% CI) [†]	0.62 (0.43–0.88)	

	TTFields + TMZ (n = 209)	TMZ Alone (n = 95)
Median OS, months	16.9	14.7
Range, months	9.7–28.2	9.8–24.8
HR (95% CI) [†]	0.66 (0.49–0.85)	



Optune has a strong safety profile with no significant increase in serious AEs compared with TMZ alone

FOR MORE INFORMATION, USE THE QR CODES:




Incidence of grade 3/4 AEs occurring in ≥5% of patients during 5 years of follow-up	Optune + TMZ (n=456) %	TMZ alone (n=216) %
≥1 AE	48	44
Blood and lymphatic system disorders Thrombocytopenia	13 9	11 5
Gastrointestinal disorders	5	4
Asthenia, fatigue, and gait disturbance	9	6
Infections	7	5
Injury, poisoning, and procedural complications (falls and medical device site reaction)	5	3
Metabolism and nutrition disorders (anorexia, dehydration, and hyperglycemia)	4	5
Musculoskeletal and connective tissue disorders	5	4
Nervous system disorders Seizures	24 6	20 6
Respiratory, thoracic, and mediastinal disorders (pulmonary embolism, dyspnea, and aspiration pneumonia)	5	5




both HCPs and patients reported stable quality of life up to 1 year of Optune use

FOR MORE INFORMATION, USE THE QR CODES:



real-world evidence showed ndGBM median overall survival extension by over 12 months in the high use TTFields group

FOR MORE INFORMATION, USE THE QR CODE:



Neuro-Oncology Advances

476 | 1-6 | 2022 | <https://doi.org/10.1093/advances/nwaa019> | Advance Access date 19 September 2022

Determinants of tumor treating field usage in patients with primary glioblastoma: A single institutional experience

Matthew T. Ballo, Kathleen W. Qualls, L. Madison Michael, Jeffrey M. Sonenson, Brandon Baughman, Saradeen Karri-Wolkoff, and Manjari Pandey

Department of Radiation Oncology, West Cancer Center & Research Institute, Memphis, Tennessee, USA (M.T.B., K.W.Q.); Neurosurgery, Sitemp Memphis Neurosurgery Clinic, Memphis, Tennessee, USA (L.M.M., J.M.S., B.B.); Department of Medical Oncology, West Cancer Center and Research Institute, Memphis, Tennessee, USA (S.K.W., M.P.)

Corresponding Author: Matthew T. Ballo, MD, Department of Radiation Oncology, West Cancer Center and Research Institute, 7985 Wolf River Blvd, Germantown, TN 38138, USA (mailto:mballo@westcancer.com)

Abstract

Background. Determinants of tumor treating fields (TTFields) usage in patients receiving combined modality therapy for primary CNS wild-type glioblastoma are currently unclear.

Methods. Ninety-one patients underwent maximal debulking surgical resection, completed external beam radiotherapy with concurrent temozolomide (TMZ), and initiated adjuvant TMZ with or without TTFields. We performed a retrospective analysis of patient, tumor, and treatment-related factors that affected TTFields usage.

Results. We identified three TTFields usage subgroups: 32 patients that declined TTFields, 40 patients that started, but had monthly compliance of less than 75% or used it for less than 2 months, and 19 patients who used TTFields for 2 or more months and maintained average monthly compliance greater than 75%. With 26.5 months median follow-up for surviving patients, the 1- and 3-year actuarial overall survival for all patients was 80% and 19%, respectively. On multivariate analysis TTFields use ($P = .03$), extent of surgical resection ($P = 0.02$), and MGMT methylation status ($P = .07$) were significantly associated with overall survival. TTFields usage was explored as a continuous variable and higher average usage was associated with longer overall survival ($P = .03$). There was no relationship between patient, tumor, or treatment-related factors and a patient's decision to use TTFields.

Conclusions. No subgroup of patients was more or less likely to initiate TTFields therapy and no subgroup was more or less likely to use TTFields as prescribed. The degree of TTFields compliance may be associated with improved survival independent of other factors.

Key Point


- It is reasonable to offer all patients with primary glioblastoma TTFields therapy as we could not identify a group that was more or less likely to discontinue therapy or unable to initiate therapy. Patients benefit from TTFields regardless of tumor or patient characteristics.

Glioblastoma is the most common and aggressive primary malignant brain tumor diagnosed in adults and has a poor survival has historically been only 14.6 months.^{1,2} Tumor treating fields (TTFields) represent a novel therapy in the treatment of glioblastoma. TTFields deliver low-intensity, intermediate-frequency (200 kHz) alternating electric fields and temozolomide (TMZ) chemotherapy, median overall survival has historically been only 14.6 months.^{1,2} Tumor treating fields (TTFields) represent a novel therapy in the treatment of glioblastoma. TTFields deliver low-intensity, intermediate-frequency (200 kHz) alternating electric fields

© The Author(s) 2022. Published by Oxford University Press, the Society for Neuro-Oncology and the European Association of Neuro-Oncology. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted reuse, distribution, and reproduction in any medium, provided the original work is properly cited.

real-world evidence validates EF-14 with statistically significant improvement in PFS and OS in Chinese patients with ndGBM

FOR MORE INFORMATION, USE THE QR CODE:



	Median OS	Median PFS
TTFields + TMZ	21.8 mo	16.0 mo
TMZ alone	15.0 mo	11.0 mo

Journal of Clinical Medicine MDPI

Article
Tumor Treating Fields Combine with Temozolomide for Newly Diagnosed Glioblastoma: A Retrospective Analysis of Chinese Patients in a Single Center

Chen et al. 2024, Volume 11, 5855. <https://doi.org/10.3390/jcm11195855>

Abstract: Introduction: TTFields plus temozolomide (TTFields/TMZ) extended survival versus TMZ alone in newly diagnosed glioblastoma (GBM) patients in the EF-14 trial. We have reported a retrospective analysis of newly diagnosed Chinese GBM patients who received TTFields/TMZ treatment and TMZ treatment from August 2019 to May 2023 in Huashan Hospital in Shanghai, Shanghai, China. Overall survival (OS) and progression-free survival (PFS) rates were constructed using the Kaplan-Meier method. A Cox proportional hazards regression model, propensity score matched data, and inverse probability of treatment weighting (IPTW) based on propensity score were used to assess the effect of TTFields and account for confounding factors. Results: In the preliminary analysis, the median PFS in the TTFields/TMZ group was 16 months (95% CI: 9.4–24.0) versus 11 months (95% CI: 6.52 to TMZ group $p < 0.05$). Median overall survival was 21.8 months (95% CI: 17.4–26.1) with TTFields/TMZ versus 15 months (95% CI: 11.48 to TMZ alone. The multivariate analysis identified age, sex, KPS score, EBV status, and TTFields use as clinically prognostic factors. After PFS adjustment, the results among the groups were similar, except that the multivariate use of MGMT promoter remained high in the TMZ group (12/17 months $p < 0.01$). Upon IPTW survival analysis, TTFields was associated with a significantly lower risk of death (HR = 0.39 in OS, 95% CI: 0.20–0.67) and progression (HR = 0.20, 95% CI: 0.14–0.29) compared with TMZ group. Conclusion: In the final analysis of our single-center Chinese patients with glioblastoma, adding TTFields to temozolomide chemotherapy resulted in statistically significant improvement in PFS and OS.


Keywords: glioblastoma; tumor treating fields; chemotherapy; retrospective cohort

Copyright: © 2024 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<http://creativecommons.org/licenses/by/4.0/>).

J. Clin. Med. **2024**, *11*, 5855. <https://doi.org/10.3390/jcm11195855> <https://www.mdpi.com/journal/jcm>

post-approval study supports safety and efficacy profile of TTFields in ndGBM Japanese patients, validating EF-14 improved survival rates

FOR MORE INFORMATION, USE THE QR CODE:



	1-year survival	2-year survival
TMZ alone	65%	31%
TTFields + TMZ	77.9%	53.6%
TTFields + TMZ high use group	89.5%	65.8%

JJCO *Japanese Journal of Clinical Oncology* 2023, 14
 https://doi.org/10.1093/jjco/hyad001
 Original Article

Safety and efficacy of tumour-treating fields (TTFields) therapy for newly diagnosed glioblastoma in Japanese patients using the Novo-TTF System: a prospective post-approval study

Ryo Nishikawa^{1,2}, Fumiyuki Yamasaki³, Yoshiki Arakawa⁴, Yoshihiro Muraguchi⁵, Yoshitaka Naita⁶, Shota Tanaka⁷, Shigeru Yamaguchi⁸, Akitsake Mukasa⁹ and Masayuki Kanomori¹⁰

¹Department of Neuro-Oncology/Neurosurgery, Satama Medical University International Medical Center, Satama, Japan, ²Department of Neurosurgery, Hiroshima University Hospital, Hiroshima, Japan, ³Department of Neurosurgery, Kyoto University Graduate School of Medicine, Kyoto, Japan, ⁴Department of Neurosurgery, Tokyo Women's Medical University Hospital, Tokyo, Japan, ⁵Department of Neurosurgery and Neuro-Oncology, National Cancer Center Hospital, Tokyo, Japan, ⁶Department of Neurosurgery, The University of Tokyo Hospital, Tokyo, Japan, ⁷Department of Neurosurgery, Hokkaido University Hospital, Sapporo, Japan, ⁸Department of Neurosurgery, Kumamoto University Hospital, Kumamoto, Japan and ⁹Department of Neurosurgery, Shinshu University Hospital, Sendai, Japan

*For reprints and all correspondence: Ryo Nishikawa, Satama Medical University International Medical Center, Department of Neuro-Oncology/Neurosurgery, 1300-1 Satama, Hokkaido-shi, Satama 069-1388, Japan. E-mail: rnohka2020@gmail.com


Received 16 July 2022; Revised 17 November 2022; Editorial Decision 28 December 2022; Accepted 2 January 2023

Abstract
Background: Tumour-treating fields therapy is a locoregional, anti-cancer treatment. Efficacy and safety of tumour-treating fields therapy in adults with newly diagnosed glioblastoma were demonstrated in the pivotal phase 3 EF-14 study (NCT02102010). Here, we report post-approval data of tumour-treating fields therapy in Japanese patients with newly diagnosed glioblastoma. **Methods:** Unlicensed post-marketing surveillance data from Japanese patients with newly diagnosed glioblastoma treated with tumour-treating fields therapy (December 2016–June 2020) were retrospectively analysed. The primary endpoints were skin, neurological and psychiatric adverse events. The secondary endpoints were 1- and 2-year overall survival rates, and the 6-month progression-free survival. Adverse events were analysed using MedDRA v24.0. The overall survival and progression-free survival were assessed using the Kaplan–Meier survival analysis (log-rank testing). The Cox proportional hazard regression analyses were also performed. **Results:** Forty patients with newly diagnosed glioblastoma were enrolled (85.2% male; median age 59 years; median baseline Karnofsky Performance Scale score 90). The most common tumour-treating fields-therapy-related adverse event was beneath-army local skin reaction (10% of patients). The adverse events were mostly mild to moderate in severity. Neurological disorders were observed in 2.5% patients (one patient reported dysesthesia). No psychiatric disorders were reported. The 1- and 2-year overall survival rates were 77.9% (95% CI: 65.1–88.3) and 53.6% (65.1–68.7%), respectively. The 6-month progression-free survival was 72.5% (95.1–87.8%). These survival

© The Author(s) 2023. Published by Oxford University Press on behalf of Japanese Society for Clinical Oncology. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted re-use, distribution, and reproduction in any medium, provided the original work is properly cited.

the most extensive study of ndGBM patients covering 18 year period confirms TTFields' positive effect on PFS and OS

FOR MORE INFORMATION, USE THE QR CODE:




	median PFS	median OS
TTFields + TMZ	19.75 mo	31.67 mo
TMZ alone (EF-14)	12.45 mo	24.80 mo




meta-analysis in ndGBM showed significant improvement in OS, and usage $\geq 75\%$ consistently prolonged survival, corroborating pivotal trial data

FOR MORE INFORMATION, USE THE QR CODE:



TTFields therapy provide consistent activity for patients with GBM irrespective of molecular alterations

FOR MORE INFORMATION, USE THE QR CODE:



Neuro-Oncology Advances

Volume 4, Issue 1, January-December 2022 | <https://doi.org/10.1093/advances/nwaa001> | Advance Access first published online 27 June 2022

Molecular alterations associated with improved outcome in patients with glioblastoma treated with Tumor-Treating Fields

Manoj Pandey, Joanne Xu, Sandeep Mittal, Jia Zeng, Michelle Saul, Santosh Kesari, Amir Azadi, Herbert Newton, Karina Deniz, Katherine Ledner, Ashley Sumral, W. Michael Kim, and Emil Lou*

West Cancer Center and Research Institute, Memphis, Tennessee, USA (M.P.); Carol Life Sciences, Phoenix, Arizona, USA (J.Z.); J.Z., M.S., W.M.K.'s Virginia Tech Carilion School of Medicine, Roanoke, Virginia, USA (M.P.); Pacific Neuroscience Institute, East Joan's Cancer Institute, Santa Monica, California, USA (S.K.); Arizona Oncology, Billings, Phoenix, Arizona, USA (K.D.); Neuro-Oncology Center, Aaker Health Cancer Institute, Chicago, Illinois, USA (H.N.); Division of Hematology, Oncology and Transplantation, Masonic Cancer Center, University of Minnesota, Minneapolis, Minnesota, USA (H.O.); E.L.'s Leander Cancer Institute, Charlotte, North Carolina, USA (E.L.)

*Corresponding Author: Emil Lou, MD, PhD, FNCP, Associate Professor of Medicine, Division of Hematology, Oncology and Transplantation, University of Minnesota, Mayo West Clinic 401, 409 Delaware Street SE, Minneapolis, MN 55455, USA (emil.lou@mayo.edu)

Abstract

Background. The genomic and overall biology landscape of glioblastoma (GB) has become clearer over the past 2 decades, as predictive and prognostic biomarkers of both de novo and transformed forms of GB have been identified. The oral chemotherapeutic agent temozolomide (TMZ) has been integral to standard-of-care treatment for nearly 2 decades. More recently, the use of non-pharmacologic interventions, such as application of alternating electric fields, called Tumor-Treating Fields (TTFields), has emerged as a complementary treatment option that increases overall survival (OS) in patients with newly diagnosed GB. The genomic factors associated with improved or lack of response to TTFields are unknown.

Methods. We performed comprehensive genomic analysis of GB tumors resected from 58 patients who went on to receive treatment using TTFields, and compared results to 57 patients who received standard treatment without TTFields.

Results. We found that molecular driver alterations in NF1, and wild-type PIK3CA and epidermal growth factor receptor (EGFR) were associated with increased benefit from TTFields as measured by progression-free survival (PFS) and OS. There were no differences when stratified by IP2 status. Meanwhile, NF1, PIK3CA, and EGFR status were combined as a Molecular Survival Score, the combination of the 3 factors significantly correlated with improved OS and PFS in TTFields treated patients compared to patients not treated with TTFields.

Conclusions. These results shed light on potential driver and passenger mutations in GB that can be validated as predictive biomarkers of response to TTFields treatment, and provide an objective and testable genomic-based approach to assessing response.

Key Points

- Alterations in NF1 were associated with increased benefit from TTFields.
- Wild-type PIK3CA and EGFR also aligned with increased benefit from this approach.
- These results provide insight into molecular differences that can be validated to tailor treatment.

© The Author(s) 2022. Published by Oxford University Press, the Society for Neuro-Oncology and the European Association of Neuro-Oncology. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted reuse, distribution, and reproduction in any medium, provided the original work is properly cited.

review article identifies TTFields therapy as one of few factors driving increased overall survival in GBM patients since the 2005 Stupp-protocol

FOR MORE INFORMATION, USE THE QR CODE:




Neff et al.

Characteristic	N	HR ¹	95% CI ¹	p-value
Age (years)	19,414	1.02	1.02, 1.03	<0.001
Sex				
Female	8,046	—	—	reference
Male	11,368	1.10	1.07, 1.14	<0.001
Elkhauser Comorbidity Score	19,414	1.01	1.01, 1.01	<0.001
Tumor-Treating Fields (ever)				
No	16,353	—	—	reference
Yes	3,061	0.77	0.73, 0.80	<0.001
Received radiation or radiosurgery (ever)				
No	7,370	—	—	reference
Yes	12,044	0.88	0.85, 0.91	<0.001
Bevacizumab (ever)				
No	15,741	—	—	reference
Yes	3,673	0.85	0.82, 0.88	<0.001

In this commercially insured dataset, TTFields improved OS to a greater extent (HR=0.77) vs. Bevacizumab (HR=0.85) or Radiation use (HR=0.88)


TTFields subset n=3,061 over 6 years

Mrugala et al.




AEs were consistent with the safety profile from the pivotal EF-11 and EF-14 clinical studies

n=23,822 over 11 years




the therapeutic potential of TFields becoming a research “hotspot”

FOR MORE INFORMATION, USE THE QR CODE:



Number of annual publications, annual cumulative number of publications and annual total citations of TFields related literature from 2007 to September 2022. (Decline in 2022 citations due to partial year)

28.5%
AVERAGE INCREASE
IN THE CUMULATIVE
NUMBER OF
PUBLICATIONS
RELATED TO TFIELDS



TTFields therapy activates inflammmasomes to induce adjuvant immunity in glioblastoma

FOR MORE INFORMATION, USE THE QR CODE:



The Journal of Clinical Investigation

RESEARCH ARTICLE

Tumor Treating Fields dually activate STING and AIM2 inflammmasomes to induce adjuvant immunity in glioblastoma

Dongling Chen,¹ Sen B. Lu,¹ Tarun S. Mukthimoon,¹ Andy-Alexandra Calinescu,¹ Mathew Sebastian,¹ Gan Jia,¹ Tiang Liu,¹ Ashley Chinnaiyan,¹ Maryam Rahmani,¹ and David S. Tian¹
 Tianjin New Drug and Process R&D Center for New Drug Therapy, Unit 1, 100 Department of Pharmacology and Medical Biotechnology Program, Institute of Pharmacology of Medicine, Canadian Institutes of Health Research, Toronto, ON, Canada

Tumor Treating Fields (TTFields), an approved therapy for glioblastoma (GBM) and malignant mesothelioma, employ noninvasive application of low-intensity intermediate-frequency, alternating electric fields to disrupt the mitotic spindle, leading to chromosome missegregation and apoptosis. Emerging evidence suggests that TTFields may also induce inflammation. However, the mechanism underlying this property and whether it can be leveraged therapeutically are unclear. Here, we report that TTFields induce focal disruption of the nuclear envelope, leading to cytosolic release of large mitochondrial clusters that intensely recruited and activated 2 major DNA sensors – cyclic GMP-AMP synthase (cGAS) and absent in melanoma 2 (AIM2) – and their cognate signaling molecules (STING and AIM2) to produce proinflammatory cytokines, type 1 interferons (IFN β), and TLR4-response genes in glioblastoma GBM models. TTFields-treated GBM cells induced cellular immunity in vivo and in vitro via cGAS- and TLR4- and AIM2-dependent means. Using single-cell and bulk RNA sequencing of peripheral blood mononuclear cells, we detected robust gene/T-cell activation of adaptive immunity in patients with GBM via a TLR4-based trajectory and identified a gene panel signature of TTFields effects on T-cell activation and clonal expansion. Collectively, these studies defined a therapeutic strategy using TTFields as cancer immunotherapy in GBM and potentially other solid tumors.


Introduction

Glioblastoma (GBM) is the most common and lethal brain cancer in adults and one of the least immunogenic tumors (1). Recent work has revealed striking immune dysregulation and functional impairment in patients with GBM. Broken systems of lymphocyte and energy and dysfunctional cytokine profiles among others, GBM tissues also possess profoundly immunosuppressed or cold tumor microenvironments (TME), characterized by tumor-infiltrating lymphocytes (TIL) and an abundance of regulatory cells, including myeloid-derived suppressor cells (MDSCs) and regulatory T cells (Treg) (2). The cold GBM TME suppresses levels of immune checkpoint proteins (3), and is further complicated by tumor cell-mediated genetic homogeneity (4). In addition, the blood brain barrier (BBB) prevents exposure of tumor-associated

antigens to immune cells and vice versa, severely limiting immunotherapeutic efforts (5). Overcoming these barriers presents a long-standing, multifaceted, immune-mediated tumor control. To “heat up” the cold GBM TME, recent efforts have focused on tumor cell genetic pathways with viral media, such as dendritic cell-based (DC)-based vaccination, immune checkpoint blockade, involving the cytokine milieu, or adopting CAR therapy to recruit tumor-specific cytotoxic T lymphocytes (CTLs) (6). However, it remains a challenge to bring a direct, cancer cell-intrinsic route to reversing the immunosuppressed state of the GBM TME.

By targeting the mitotic apparatus, and assembly of macromolecules required for the mitotic spindle structure during metaphase and the centrosome ring during anaphase, microtubules, and cytoskeleton of the cell cycle, Tumor Treating Fields (TTFields) cause chromosome missegregation and breakage and chromosome cytoplasmic separation, respectively, leading to mitotic catastrophe and cell death (7). Independent of apoptosis (8), TTFields have also been demonstrated to target the DNA damage repair and base excision (BER) (9,10) and mismatch (MMR) (11) signaling repair pathways by interfering with DNA Glnk replication (8–10) and induce chromosome structural errors during mitosis to trigger chromosome nonreciprocal cell cycle kinase-dependent aneuploidization (12). Through decreased degradation of p53 and light chain 4, p24 (LCN4) in tumor LCN4, Rb-E10. Recent reports also correlated TTFields’ ability to downregulate the plasma membrane of GBM cells, affecting particles up to 20 kDa in size per

Conflict of interest: The authors have nothing to disclose.
Copyright: © 2024 Chen et al. This is an open-access article distributed under the terms of the [Creative Commons Attribution 4.0 International license](https://creativecommons.org/licenses/by/4.0/).
Received: March 14, 2024; **Accepted:** February 16, 2024; **Published:** April 16, 2024.
Address correspondence: Dr. Tian to David.S.Tian@utoronto.ca.




there is early evidence of efficacy in newly diagnosed GBM patients when TTFields therapy is added to immune checkpoint inhibitors

FOR MORE INFORMATION, USE THE QR CODE:



Overall Survival 2-THE-TOP single arm study vs. external controls



WFNOS 2022 Top 10 Session 2 / March 26 (Sat), 10:15-11:30

BTRT

Phase 2
Phase 2 study of pembrolizumab plus TTFields plus temozolomide in patients with newly diagnosed glioblastoma (2-THE-TOP)

David Tran, Ashley Ghazieddin, Dongjiang Chen, Maryam Rahmani
Department of Neurosurgery, Division of Neuro-Oncology, University of Florida, United States

Background: Emerging data indicate that TTFields, the new anti-mitotic treatment for GBM, stimulate immunity via the type-1 interferon (IFN1) pathway of STING and AIM2 inflammasomes. Thus, we hypothesize that TTFields synergize with immune checkpoint inhibitors to induce anti-tumor immunity in GBM.

Methods: We conducted a phase 2 study combining pembrolizumab, TTFields and maintenance TMZ in 26 patients with newly diagnosed GBM (nGBM). To delineate immune effects of TTFields from pembrolizumab, TTFields was started at cycle 1 of TMZ and pembrolizumab (200 mg Q3 weeks) at cycle 2. The primary endpoint was PFS vs. the historical control of TTFields plus TMZ (JAMA/318/2006-2106) and immune signatures of TTFields and pembrolizumab by single-cell genomics of PBMCs. Secondary endpoints included toxicity and OS.


Results: As of 09/24/2022, 26 patients with a median age of 60 years were enrolled. Fourteen (54%) had biopsy only or subtotal resection. Nineteen (73%) had immunohistochemical MGMT and 3 (12%) had an IDH mutation. The median follow-up was 10 and 18.2 months for PFS and OS, respectively. Thirteen (50%) were progression free and 10 (38%) were alive. Of 22 patients with follow-up of 6 months, the median PFS was 3.11 vs. 6.7 months in the control. Six (27%) patients with measurable tumors have achieved partial to complete objective response. We assessed 65/70 PBMCs in 12 patients before and after TTFields and detected robust post-TTFields T cell activation in 31 of 12 patients via the IFN1 regulatory which was strongly correlated with TCRβ-αβ clonal expansion (Spearman coefficient $r = 0.8$, $p = 0.014$). Importantly, we defined a T cell-based gene signature of TTFields effect on TCRβ clonal expansion. The most common network adverse events were thrombocytopenia, asthenia, and metabolic disturbances in 4 (15%), 3 (12%), and 2 (8%) patients, respectively.

Conclusions: The triple combination is well tolerated and shows early evidence of efficacy in nGBM patients. Survival and molecular data will be updated.

Keywords: TTFields immunotherapy; pembrolizumab; STING; single cell analysis

STELLAR phase 2 trial evaluated TTFields therapy + pemetrexed and cisplatin or carboplatin in MPM

N = 80



Start date: February 2015
Primary completion: April 2018
Study completion: April 2018
Study sites: 13 (Europe)

Primary endpoints:


- OS

Secondary endpoints:


- PFS, ORR (modified RECIST criteria for MPM), safety

MPM patients who used Optune Lua first line achieved 18.2 months median overall survival


FOR MORE INFORMATION, USE THE QR CODE:



Adapted from Ceresoli GL et al. 2019



EF-15 phase 2 trial evaluated TTFields therapy + pemetrexed in NSCLC



Start date: May 2008

Primary completion: July 2011

Study completion: July 2011

Study sites: 4 (Switzerland)

Primary endpoints:


- Device related toxicity (P1), Time to in-field progression (P2)

Secondary endpoints:

- OS, ORR, time to systemic progression, safety

TTFields therapy together with pemetrexed improved disease control within the treatment field in second line NSCLC


FOR MORE INFORMATION, USE THE QR CODE:



	Median in-field PFS	Median PFS	Median OS	1yr Survival
TTFields + Pemetrexed	6.5 mo	5.0 mo	13.8 mo	57.0%
Pemetrexed alone	n/a	2.9 mo	8.3 mo	29.7%



LUNAR phase 3 trial evaluated TTFields therapy + SOC in metastatic NSCLC, post-platinum



Start date: December 2016

Primary completion: December 2022

Study completion: December 2022

Study sites: 124

Primary endpoints:

- OS

Secondary endpoints:


- OS (by cohort), PFS, ORR, QoL, safety

TTFields therapy together with either standard of care therapies or immune checkpoint inhibitor improved overall survival in second-line NSCLC



FOR MORE INFORMATION, USE THE QR CODE:




Overall survival (ITT population)



Overall survival (ICI-treated patients)



PANOVA phase 2 trial evaluated TTFields therapy + gemcitabine +/- nab-paclitaxel in pancreatic cancer



Start date: Nov 2013

Primary completion date: Dec 2017

Study completion date: Dec 2017

Study sites: 6 (Europe)

Primary endpoint:


- Safety

Secondary endpoints:

- TTFields monthly usage, PFS, OS


TTFields therapy together with chemotherapy were well tolerated for patients with advanced pancreatic cancer

FOR MORE INFORMATION, USE THE QR CODE:



OS	
Median, mo	NR
95% CI	8.4–NA
1-year survival	72%

	Median PFS	Median OS	One-year Survival	Partial Response Rate	Stable Disease
TTFields + gemcitabine	8.3 mo	14.9 mo	55%	30%	30%
gemcitabine alone	3.7 mo	6.7 mo	22%	7%	28%
TTFields + gemcitabine + nab-paclitaxel	12.7 mo	Not yet reached	72%	40%	47%
gemcitabine + nab-paclitaxel alone	5.5 mo	8.5 mo	35%	23%	27%




encouraging response rate and durability signals in EF-31 phase 2 gastric cancer trial

FOR MORE
INFORMATION, USE
THE QR CODE:



EF-31 PHASE 2 PILOT TRIAL DESIGN¹




	OBJECTIVE RESPONSE RATE	MEDIAN PROGRESSION-FREE SURVIVAL	DURATION OF RESPONSE	ONE-YEAR SURVIVAL
TTFields + chemotherapy	50%	7.8mo	10.3mo	72%
SOC chemotherapy ²	41-45%	6.9mo	6.9mo	48%

encouraging signals in liver cancer despite poor prognosis and low treatment exposure in HEPANOVA phase 2 trial

FOR MORE INFORMATION, USE THE QR CODE:



HEPANOVA PHASE 2 PILOT TRIAL DESIGN²



76%

DISEASE CONTROL RATE (n=21)

VS. 43% CONTROL³

9.5%

OBJECTIVE RESPONSE RATE (n=21)

VS. 4.5% CONTROL

91%


DISEASE CONTROL RATE

18%

OBJECTIVE RESPONSE RATE

patients that received ≥ 12 wks of TTFields (n=11)

platform technology driving robust clinical pipeline






tumor treating fields mechanism of action overview

patients with aggressive solid tumors often face suboptimal survival outcomes, despite advancements in treatment modalities

These outcomes are due to diverse treatment challenges, including:



Therapeutic tumor resistance



Drug-to-drug interactions




Additive systemic toxicities

With a poor survival outlook, physicians and patients need additional treatment strategies

Tumor Treating Fields (TTFields) are electric fields that exert physical forces to kill cancer cells via a variety of mechanisms



TTFields spare healthy cells because they have different properties than cancer cells across a range of tumor types



a growing body of evidence supporting multiple mechanisms of action

FOR MORE
INFORMATION, USE
THE QR CODE:



- Preclinical research has shown interference with cancer cell motility and migration, activation of anti-tumor immunity, downregulation of genes important for DNA damage repair, and other potential mechanisms
- May demonstrate enhanced effects across solid tumor types when used with chemotherapy, radiotherapy, immune checkpoint inhibition, or PARP inhibition in preclinical models

JOURNAL ARTICLE ACCEPTED MANUSCRIPT

Anti-cancer mechanisms of action of therapeutic alternating electric fields (tumor treating fields [TTFields])

Shadi Shams, Chirag B Patel


Journal of Molecular Cell Biology, mjac047, <https://doi.org/10.1093/jmcb/mjac047>
Published: 16 August 2022

PDF Split View Cite Permissions Share

Abstract



Despite improved survival outcomes across many cancer types, the prognosis remains grim for certain solid organ cancers including glioblastoma and pancreatic cancer. Invariably in these cancers, the control achieved by time-limited interventions such as traditional surgical resection, radiation therapy, and chemotherapy is short-lived. A new form of anti-cancer therapy called therapeutic alternating electric fields (AEFs) or tumor treating fields (TTFields) has been shown, either by itself or in combination with chemotherapy, to have anti-cancer effects that translate to improved survival outcomes in patients. Although the pre-clinical and clinical data are promising, the mechanisms of TTFields are not fully elucidated. Many investigations are underway to better understand how and why TTFields is able to selectively kill cancer cells and impede their proliferation. The purpose of this review is to summarize and discuss the reported mechanisms of action of TTFields from pre-clinical studies (both *in vitro* and *in vivo*). An improved understanding of how TTFields works will guide strategies focused on the timing and combination of TTFields with other therapies, to further improve survival outcomes in patients with solid organ cancers.

Tumor Treating Fields have multiple, distinct mechanisms of action




TTFields have been shown to disrupt mitosis in cancer cells by exerting physical forces on their polar components

FOR MORE INFORMATION, USE THE QR CODE:



TTFields have been shown to alter the organization and dynamics of the cytoskeleton, disrupting cancer cell motility and migration

FOR MORE INFORMATION, USE THE QR CODE:



a) microtubule; b) TTFields; c) tubulin aligned with field; j) actin fiber; k) integrin; l) focal adhesion; m) extracellular matrix.

A model illustrating the mechanism by which TTFields modulates cancer cell motility.

- (1) Microtubules are required to specify the direction of cell movement. GEF-H1 catalytic activity is downregulated through microtubule binding.
- (2) TTFields exert directional forces on polar tubulins leading to their alignment in the direction of the field. This, in turn, leads to the reorganization of the microtubule network resulting in changes in the abundance of microtubules and initiation of the GEF-H1/RhoA/ROCK signaling pathway
- (3) to increase actin bundling
- (4) and formation of focal adhesions,
- (5) which disrupt cell polarity and migration directionality.

cancers
Tumor Treating Fields (TTFields) Hinder Cancer Cell Motility through Regulation of Microtubule and Actin Dynamics


Tali Voloshin¹, Roni Sass Schneiderman¹, Alexandra Volodin, Reuben Eady Shavit, Noa Karzen, Erez Zevi, Ilan Keren, Anat Kater-Goldberg, Roni Yan, Meha Gilad¹, Zeev Ronson, Uri Weising and Nissim Katz¹
 November 15th, Tumor Treating Fields (TTFields) application modulates cancer cell motility through regulation of microtubule and actin dynamics. *Cancers* 2020;12(10):1-18. doi:10.3390/cancers12103016

Simple Summary: Tumor Treating Fields (TTFields), encompassing alternating electric fields within the intermediate frequency range, is an anticancer treatment delivered to the tumor region through transducer arrays placed non-invasively on the skin. Although established as an anti-neurotic treatment modality, the anti-metastatic potential of TTFields and their effect on rapid cytoskeletal dynamics during cellular motility warrant further investigation. In this study, we report that TTFields application induces changes in microtubule organization leading to interference with the directionality and abundance of cancer cell migration. We show that these changes in microtubule organization result in activation of GEF-H1/RhoA/ROCK signaling pathway, and the consequent formation of focal adhesions and changes in actin cytoskeleton architecture. Together, these results propose a novel mechanism by which TTFields induce changes in microtubule and actin organization and dynamics, thereby disrupting processes important for polarity generation and motility in cancer cells.

Abstract: Tumor Treating Fields (TTFields) are non-invasive, alternating electric fields within the intermediate frequency range (100–300 kHz) that are utilized as an anticancer cancer treatment. TTFields are transcutaneously delivered to the tumor region through 2 pairs of transducer arrays placed on the skin. This novel treatment modality has been FDA-approved for use in patients with glioblastoma and malignant pleural mesothelioma but has not been clinically tested in metastatic efficacy and safety, and is currently under investigation in other types of solid tumors. TTFields were shown to induce an anti-mitotic effect by exerting bidirectional forces on highly polar intracellular organelles, such as tubulin and septin microtubules, causing abnormal microtubule polymerization during spindle formation as well as aberrant cleavage furrow formation. Previous studies have demonstrated that TTFields inhibit metastatic progression in cancer cells. However, the consequences of TTFields application on cytoskeleton dynamics remain undetermined. In this study, methods utilized in combination to study the effects of TTFields on cancer cell motility through regulation of microtubule and actin dynamics included confocal microscopy, computational tools, and biochemical analyses. Mechanisms by which TTFields treatment disrupted cellular polarity were (1) interference with microtubule assembly and directionality; (2) altered regulation of G-actin microtubule exchange factor (GEF-H1), Rho GTPase family member A (RhoA), and Rho-associated coiled-coil kinase (ROCK) activity; and (3) induced formation of radial protrusions of peripheral actin filaments and focal adhesions. Overall, these data identified discrete effects of TTFields that disrupt processes crucial for cancer cell motility.

TTFields-mediated cell disruption activates the immune system and triggers a downstream antitumor cell response

FOR MORE INFORMATION, USE THE QR CODE:



TTFields induces downstream immunogenic cell death, including release of DAMPs (damage-associated molecular patterns)

Cancer Immunology, Immunotherapy (2020) 69(7):1191–1204
<https://doi.org/10.1007/s00262-020-02034-7>

ORIGINAL ARTICLE

Tumor-treating fields (TTFields) induce immunogenic cell death resulting in enhanced antitumor efficacy when combined with anti-PD-1 therapy

Tali Velichko¹, Nava Kaynan¹, Shai Davidi¹, Naara Porel¹, Anna Steingass¹, Ronit S. Schneiderman¹, Elmar Zeevi¹, Mital Munster¹, Roni Blat¹, Catherine Tempel Breen¹, Shay Cahal¹, Arnon Itzhaki¹, Moshé Glid¹, Eliran G. Krizan¹, Avi Mendelsohn¹, Adirah Khouf¹, Yoram Kafri¹

Received: 3 September 2019 / Accepted: 24 February 2020 / Published online: 6 March 2020
 © The Author(s) 2020

Abstract
 Tumor-treating fields (TTFields) are alternating electric fields in a specific frequency range (100–300 kHz) delivered to the human body through transducer arrays. In this study, we evaluated whether TTFields-mediated cell death can elicit antitumoral immunity and hence would be effectively combined with anti-PD-1 therapy. We demonstrate that in TTFields-treated cancer cells, damage-associated molecular patterns including high-mobility group B1 and adenosine triphosphatase are released and colocalize in a pool on the cell surface. Moreover, we show that TTFields treatment promotes the engagement of cancer cells by dendritic cells (DCs) and DCs maturation in vivo, as well as recruitment of immune cells in vivo. Additionally, our study demonstrates that the combination of TTFields with anti-PD-1 therapy results in a significant decline

This work was presented at a poster at the annual meeting of the American Association of Immunologists (AAI), May 13–17, 2016, Seattle, WA, USA (P0797); the annual meeting of the American Association for Cancer Research (AACR), May 1–5, 2017, Washington, DC, USA (100); the annual meeting of the Society for Neuro-Oncology (SNO), November 10–15, 2017, San Francisco, CA, USA (108); the annual meeting of the Society for Neuro-Oncology (SNO), November 10–15, 2018, New Orleans, LA, USA (108); the European Association of Neuro-Oncology (EANO) Meeting, October 10–12, 2018, Stockholm, Sweden (P0742); the annual meeting of the American Association for Cancer Research (AACR), March 26–April 1, 2019, Atlanta, GA, USA (108); the International Anti-Cancer Immunotherapy Conference (CACI), November 21–26, 2019, Paris, France (108); the annual meeting of the Society for Neuro-Oncology (SNO), November 10–15, 2019, National Harbor, MD, USA (108); the annual meeting of the American Society for Radiation Oncology (ASTRO), September 13–16, 2019, Chicago, IL, USA (108); the Breast Cancer Meeting, May 21–24, 2019, Seattle, Washington, WA (108); the Lung Cancer Congress (LCC), April 10–13, 2019, Geneva, Switzerland (108); the Multidisciplinary Thoracic Cancer Symposium, March 14–16, 2019, San Diego, CA, USA (108).

Tali Velichko and Nava Kaynan have contributed equally to this work.

Electronic supplementary material The online version of this article (<https://doi.org/10.1007/s00262-020-02034-7>) contains supplementary material, which is available to authorized users.


Extended author information available on the last page of the article.

TTFields downregulate genes important for DNA damage repair

FOR MORE INFORMATION, USE THE QR CODE:



- TTFields disrupt DNA damage repair in cancer cells by downregulating genes that are part of the well-known FA-BRCA pathway^{1,2}



OPEN

Tumor-treating fields elicit a conditional vulnerability to ionizing radiation via the downregulation of BRCA1 signaling and reduced DNA double-strand break repair capacity in non-small cell lung cancer cell lines

Narasimha Kumar Karanam¹, Kalyanram Srinivasan¹, Lingshu Ding¹, Brock Blain¹, Debabrata Saha¹ and Michael D Story^{1*}

Cell Death and Disease (2017) 8:4271 | doi:10.1038/s41419-017-3136-1

The use of tumor-treating fields (TTFields) has revolutionized the treatment of recurrent and newly diagnosed glioblastoma (GBM). TTFields are low-frequency, intermediate-frequency alternating electric fields that are applied to tumor regions and cells using non-invasive arrays. The predominant mechanism by which TTFields are thought to kill tumor cells is the disruption of mitosis. Using the non-small cell lung cancer (NSCLC) cell lines we found that there is a variable response to cell proliferation and cell killing between the NSCLC cell lines that are independent of cell status. TTFields treatment increased the G2M population, with a concomitant reduction in S-phase cells followed by the appearance of a sub-G1 population indicative of apoptosis. Temporal changes in gene expression during TTFields exposure were evaluated to identify molecular signaling changes underlying the differential TTFields response. The most differentially expressed genes were associated with the cell cycle and cell proliferation pathways. However, the expression of genes found within the BRCA1-DNA damage response were significantly downregulated (P < 0.05) during TTFields treatment. DNA double-strand break (DSB) repair foci increased when cells were exposed to TTFields as did the appearance of chromosome aberrations, suggesting an interference mechanism responsible for cell killing involving DNA repair. Exposing cells to TTFields immediately following ionizing radiation resulted in increased chromosomal aberrations and a reduced capacity to repair DNA DSBs, which were fully reversible in at least a portion of the enhanced cell killing seen with the combination. These findings suggest that TTFields induce a state of BRCA1 loss leading to a conditional susceptibility resulting in enhanced sensitivity to ionizing radiation and provide a strong rationale for the use of TTFields as a combined modality therapy with radiation or other DNA-damaging agents.

Cell Death and Disease (2017) 8:4271 | doi:10.1038/s41419-017-3136-1 published online 26 March 2017

Lung cancer is the second most prevalent cancer and the leading cause of cancer-related death in the United States.¹ Non-small cell lung cancer (NSCLC) is the most prevalent type, accounting for ~85% of new cases.² A plethora of treatment options, such as including surgical resection, chemotherapy, radiation therapy and immunotherapy,^{3–7} have not achieved rates for patients with stage I and II NSCLC as high as 50% and 30%, respectively. However, despite the myriad of options, 5-year survival rates for patients with late-stage NSCLC are at 10–15%.^{8,9} Thus, the need for novel treatment strategies that can be utilized alone or in combination with conventional therapies to increase survival rates.

The advent of Tumor-Treating Fields (TTFields), a novel physical treatment modality, has been effective for the treatment of adult, recurrent glioblastoma (GBM) and recurrent tumors.^{10–12} TTFields are non-invasive and deliver a non-thermal, non-ionizing, intermediate frequency (100–200 kHz) alternating electric field across the tumor base.¹³ TTFields create a heterogeneous interfacial environment that induces a dielectrophoretic movement of polar molecules towards the region of highest field intensity, effectively generating asymmetric and other critical biochemical functions.¹⁴ As such, TTFields perturb and target cancer cells through the inhibition of cell proliferation, effectively opening non-dividing non-mitotic cells to cell death. TTFields do not ablate nerves and muscle because of their high frequency, and do not penetrate the brain due to their non-invasive nature. The use of TTFields in combination with chemotherapy, immunotherapy, and/or radiation therapy, for the treatment of recurrent and newly diagnosed glioblastoma (GBM) is currently under investigation.¹⁵ Clinical trials are ongoing or recruiting for cancer and glioblastoma (www.novocure.com).

TTFields are known to disrupt cellular proliferation and induce apoptotic pathways in dividing cancer cells across a variety of human solid tumor types (cell lines).¹⁶ Prevention of proper formation of the mitotic spindle apparatus and the inhibition of other mitotic structures have been proposed as the mechanism by which TTFields kill dividing cells.^{17–19} Specifically, TTFields require stable microtubule depolymerization and the miscalculation of length. This results in

TTFields is a highly versatile first-in-class treatment modality

FOR MORE
INFORMATION, USE
THE QR CODE:




- TTFields therapy has significant potential for broad applicability across solid tumor types and lines of therapy
- Investigation of TTFields therapy is ongoing across clinical trials in multiple tumor types
- In approved indications, TTFields therapy is well tolerated, suggesting a low risk of additive systemic toxicity when used with other cancer treatment modalities




TTFields therapy can be added to cancer treatment modalities in approved indications


TTFields **demonstrate enhanced effects** across multiple solid tumor types, when used concomitantly with each of the following:




Chemotherapy



Radiation therapy (RT)



Targeted therapies



Immuno-oncologic (IO) agents