




tumor treating fields clinical evidence

EF-14 Phase 3 pivotal trial evaluated Optune Gio + TMZ in 695 patients with ndGBM





Start date: June 2009
Primary completion: December 2016
Study completion: March 2017
Study sites: 83 (global)

Primary endpoint:
 • Progression-free survival
Secondary endpoints:
 • Overall survival


in ndGBM, Optune Gio + TMZ provided an unprecedented long-term survival benefit

FOR MORE INFORMATION, USE THE QR CODE:



more time on Optune Gio predicted increased significant survival benefit

FOR MORE INFORMATION, USE THE QR CODE:



29.3%
vs. 4.5%

5-YEAR PROBABILITY OF SURVIVAL WITH 90% COMPLIANCE (n=43) VS SURVIVAL WITH TMZ ALONE




higher TTFields therapy dose can lead to increased efficacy

FOR MORE INFORMATION, USE THE QR CODES:




EXISTING ARRAYS AP channel, 1,364 mAmps





VS.

NEW ARRAYS AP channel, 1,685 mAmps




all analyzed subgroups experienced a benefit when adding Optune Gio to TMZ

FOR MORE INFORMATION, USE THE QR CODE:



Optune Gio was associated with increased survival in patients 65 years and older

FOR MORE INFORMATION, USE THE QR CODE:



frontiers
in Oncology

CLINICAL TRIAL
PUBLISHED IN FRONTIERS IN ONCOLOGY
DOI: 10.3389/fonc.2022.902929

Efficacy and Safety of Tumor Treating Fields (TTFields) in Elderly Patients with Newly Diagnosed Glioblastoma: Subgroup Analysis of the Phase 3 EF-14 Clinical Trial

Zul Bahri^{1*}, Chien-Feng Hsin², Andreas F. Hottinger³, Ahmed Ibrahim⁴, Gorth Nicholas⁵ and Jay-Angiang Zhu⁶

¹Department of Neurosurgery, Tri-Alice Medical Center and Tri-Alice University School of Medicine, Tri-Alice, Taiwan; ²Department of Neurosurgery, Sun Yat-sen Memorial Hospital, Sun Yat-sen University School of Medicine, Guangzhou, China; ³Department of Clinical Neurosciences, CHUQ (Université Laval) and Université de Québec, Québec, Québec, Canada; ⁴Service de Neurologie, Hôpital Ste-Justine, Université de Québec, Québec, Québec, Canada; ⁵Department of Neurology, University of Illinois at Chicago, Chicago, IL, United States; ⁶Department of Neurosurgery, University of Texas Health Science Center at Houston, Houston, TX, United States

OPEN ACCESS

Edited by: Matthew P. Froscher, Medical University of Vienna, Austria

Reviewed by: Massimo Ghisletta, Università del Piemonte Orientale, Italy; M. S. Subramanian, Cancer Research and Biotechnology Institute, India; Venkata Subrahmanyam, Venkata Institute of Technology, India

***Correspondence:** Zul Bahri, zbahri@triai.com.tw

Specialty section: This article was submitted to Neuro-Oncology and Cancer, a specialty of Frontiers in Oncology.

Received: 04 February 2022
Accepted: 27 August 2022
Published: 27 September 2022

Citation: Bahri Z, Hsin C-F, Hottinger A, Ibrahim A, Nicholas G and Zhu J-A (2022) Efficacy and Safety of Tumor Treating Fields (TTFields) in Elderly Patients with Newly Diagnosed Glioblastoma: Subgroup Analysis of the Phase 3 EF-14 Clinical Trial. *Front. Oncol.* 12:902929. doi: 10.3389/fonc.2022.902929

Keywords: glioblastoma, elderly patients, tumor treating fields, temozolomide, overall survival, safety

Background: Understudied elderly patients comprise a large segment of high-risk patients with glioblastoma (GBM) that are challenging to treat. Tumor Treating Fields (TTFields) is a locoregional, noninvasive, antimitotic therapy delivering low-intensity, intermediate-frequency, alternating electric fields to the tumor. In the phase 3 EF-14 clinical trial, TTFields (200 kHz) improved median progression-free survival (PFS) and median overall survival (OS) in patients with newly diagnosed GBM (nGBM) when added concomitantly to maintenance temozolomide (TMZ). This EF-14 subgroup analysis evaluated the safety and efficacy of TTFields in elderly patients.


Methods: All 184 patients who are ≥65 years of age were included (TTFields/TMZ combination, n=89; TMZ monotherapy, n=95; 2:1 ratio of randomization). PFS/OS were analyzed using Kaplan-Meier methodology (α=0.05). Health-related quality-of-life (HRQL) was assessed using the European Organization for Research and Treatment of Cancer (EORTC) quality-of-life questionnaire (QLQ-C15) supplemented with the brain tumor module (QLQ-BRAND). Adverse events (AEs) were evaluated using Common Terminology Criteria for Adverse Events (CTCAE) v4.0.

Results: The PFS was 6.6 months in patients randomized to the treatment group with TTFields/TMZ combination versus 3.9 months in patients treated with TMZ monotherapy (HR, 0.47; 95% CI, 0.33-0.74; P=0.0028). The OS was 17.4 months in patients treated with TTFields/TMZ combination versus 13.7 months in patients treated with TMZ monotherapy (HR, 0.51; 95% CI, 0.33-0.77; P=0.0046). Annual survival rates with TTFields/TMZ versus TMZ monotherapy were 50% (95% CI, 33-57%) versus 27% (95% CI, 18-41%; P=0.072) at 2 years, 19% (95% CI, 11-29%) versus 11% (95% CI, 4-22%; P=0.138) at 3 years, and 19% (95% CI, 7-29%) versus 0% at 5 years, respectively. There were no significant differences between groups in the preselected

Frontiers in Oncology | www.frontiersin.org September 2022 | Volume 12 | Article 902929



biopsy-only patients using Optune Gio had longer median overall survival

FOR MORE INFORMATION, USE THE QR CODE:



survival benefit occurred independently of MGMT methylation status

FOR MORE INFORMATION, USE THE QR CODE:



	TTFields + TMZ (n = 137)	TMZ Alone (n = 77)
Median OS, months	31.6	21.2
Range, months	21.1–48.5	12.3–37.9
HR (95% CI) ¹	0.62 (0.43–0.88)	

	TTFields + TMZ (n = 209)	TMZ Alone (n = 95)
Median OS, months	16.9	14.7
Range, months	9.7–28.2	9.8–24.8
HR (95% CI) ¹	0.66 (0.49–0.85)	



Optune Gio has a strong safety profile with no significant increase in serious AEs compared with TMZ alone

FOR MORE INFORMATION, USE THE QR CODES:




Incidence of grade 3/4 AEs occurring in ≥5% of patients during 5 years of follow-up	Optune + TMZ (n=456) %	TMZ alone (n=216) %
≥1 AE	48	44
Blood and lymphatic system disorders Thrombocytopenia	13 9	11 5
Gastrointestinal disorders	5	4
Asthenia, fatigue, and gait disturbance	9	6
Infections	7	5
Injury, poisoning, and procedural complications (falls and medical device site reaction)	5	3
Metabolism and nutrition disorders (anorexia, dehydration, and hyperglycemia)	4	5
Musculoskeletal and connective tissue disorders	5	4
Nervous system disorders Seizures	24 6	20 6
Respiratory, thoracic, and mediastinal disorders (pulmonary embolism, dyspnea, and aspiration pneumonia)	5	5





HCPs and patients reported stable quality of life up to 1 year of Optune Gio use

FOR MORE
INFORMATION, USE
THE QR CODES:




meta-analysis in ndGBM showed significant improvement in OS, and usage $\geq 75\%$ consistently prolonged survival, corroborating pivotal trial data

FOR MORE INFORMATION, USE THE QR CODE:



TTFields therapy provide consistent activity for patients with GBM irrespective of molecular alterations

FOR MORE INFORMATION, USE THE QR CODE:



Neuro-Oncology Advances

Volume 4, Issue 1, January-December 2022 | <https://doi.org/10.1093/advances/nwac096> | Advance Access Date 27 June 2022

Molecular alterations associated with improved outcome in patients with glioblastoma treated with Tumor-Treating Fields

Manoj Pandey, Joanne Xu, Sandeep Mittal, Ju Zeng, Michelle Saul, Santosh Kesari, Amir Azad, Herbert Newton, Karina Deniz, Katherine Ladner, Ashley Sumral, W. Michael Kim, and Emil Lau*

West Cancer Center and Research Institute, Memphis, Tennessee, USA (M.P.); Dana Life Sciences, Phoenix, Arizona, USA (J.Z.); J.Z., M.P., W.M.K.'s Virginia Rich Center School of Medicine, Roanoke, Virginia, USA (J.Z.); Pacific Neuroscience Institute, Siteman Cancer Institute, Siteman Cancer Center, California, USA (S.M.); Arizona Oncology, Billings, Phoenix, Arizona, USA (K.D.); Neuro-Oncology Center, Aaker Health Center Institute, Orlando, Florida, USA (H.N.); Division of Hematology, Oncology and Reproductive, Maimon Cancer Center, University of Minnesota, Minneapolis, Minnesota, USA (E.L.); E.L.'s Cancer Center Institute, Charlotte, North Carolina, USA (E.L.)

*Corresponding Author: Emil Lau, MD, PhD, FACP, Associate Professor of Medicine, Division of Hematology, Oncology and Reproductive, University of Minnesota, Mayo West Code 480, 430 Delaware Street SE, Minneapolis, MN 55455, USA (emil.lau@umn.edu)

Abstract

Background. The genomic and overall biology landscape of glioblastoma (GB) has become clearer over the past 2 decades, as predictive and prognostic biomarkers of both de novo and transformed forms of GB have been identified. The oral chemotherapeutic agent temozolomide (TMZ) has been integral to standard-of-care treatment for nearly 2 decades. More recently, the use of non-pharmacologic interventions, such as application of alternating electric fields, called Tumor-Treating Fields (TTFields), has emerged as a complementary treatment option that increases overall survival (OS) in patients with newly diagnosed GB. The genomic factors associated with improved or lack of response to TTFields are unknown.

Methods. We performed comprehensive genomic analysis of GB tumors resected from 58 patients who went on to receive treatment using TTFields, and compared results to 37 patients who received standard treatment without TTFields.

Results. We found that molecular driver alterations in NF1, and wild-type PIK3CA and epidermal growth factor receptor (EGFR), were associated with increased benefit from TTFields, as measured by progression-free survival (PFS) and OS. There were no differences when stratified by IP2 status. Wild-type EGFR, PIK3CA, and EGFR status were combined as a Molecular Survival Score, the combination of the 3 factors significantly correlated with improved OS and PFS in TTFields-treated patients compared to patients not treated with TTFields.

Conclusions. These results shed light on potential driver and passenger molecules in GB that can be validated as predictive biomarkers of response to TTFields treatment, and provide an objective and testable genomic-based approach to assessing response.


Key Points

- Alterations in NF1 were associated with increased benefit from TTFields.
- Wild-type PIK3CA and EGFR also aligned with increased benefit from this approach.
- These results provide insight into molecular differences that can be validated to tailor treatment.

© The Author(s) 2022. Published by Oxford University Press, the Society for Neuro-Oncology, and the European Association of Neuro-Oncology. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

the therapeutic potential of TTFields therapy becoming a research “hotspot”

FOR MORE INFORMATION, USE THE QR CODE:



Number of annual publications, annual cumulative number of publications and annual total citations of TTFields related literature from 2007 to September 2022. (Decline in 2022 citations due to partial year)

28.5%
AVERAGE INCREASE
IN THE CUMULATIVE
NUMBER OF
PUBLICATIONS
RELATED TO TTFIELDS

frontiers | Frontiers in Oncology

OPEN ACCESS

EDITED BY
Michael Card,
University of North Carolina at Chapel Hill, United States

REVIEWED BY
Yuan Zhang,
University of North Carolina at Chapel Hill, United States

*CORRESPONDENCE
Xue Du ^{1,2}, Chunbao Chen ¹, Yu Xiao ¹, Yu Cui ^{1,3}, Lu Yang ^{1,4}, Xiaochun Li ^{1,5}, Jingling Lu ^{1,6}, Ruiti Wang ^{1,7} and Bingtuan Tan ^{1,8} ^{*}

¹Department of Oncology, Affiliated Hospital of North Sichuan Medical College, North Sichuan Medical College, Nanchong, Sichuan, China; ²Department of Clinical Medicine, North Sichuan Medical College, Nanchong, Sichuan, China; ³Department of Radiotherapy, Fudan Cancer Hospital, Fudan University, Shanghai, China

RECEIVED 10 November 2022
ACCEPTED 02 November 2022
PUBLISHED 10 November 2022

Research on application of tumor treating fields in glioblastoma: A bibliometric and visual analysis

Background: Glioblastoma, one of the common tumors of the central nervous system (CNS), is prone to recurrence even after standard treatment protocols. As an innovative physiotherapy method emerging in recent years, the tumor treating fields (TTFields) technique has been approved for the treatment of glioblastoma due to its non-invasive and portable features. The purpose of this study is to evaluate and analyze the scientific trends and research trends in TTFields therapy for glioblastoma.


Methods: Publications related to TTFields therapy for glioblastoma were searched in the Web of Science Core Collection (WOSCC) database in September 2022. A bibliometric and visual analysis of publications in this field was performed mainly using CiteSpace and R software for country/region, author, journal, reference and keywords.

Results: A total of 458 publications in this field were retrieved, and 348 were finally obtained according to the search criteria, including 159 articles (64.11%) and 89 reviews (27.9%). The cumulative number of publications increased year by year, with an average growth rate (AGR) of 28.50%. The best result of Pearson correlation coefficient showed a high positive correlation between publications and citations (r=0.932, p<0.001). The USA had the largest number of publications (123, 49.80%), followed by Germany (52, 12.90%) and China (30, 12.30%). As for the country/region collaborations, the USA cooperated most closely with other countries/regions, followed by Germany and China. The degree of collaboration (DC) between countries/regions was 25.05%. The institutions with the largest number of publications were Fudan Univ [10], Harvard Med Sch [10] and Novocure [10]. Moreover, Wang F (8) possessed the greatest number of publications, followed by Wenkang L (1) and Wang F (1). The DC between authors was 37.58%. STUPP (125) was the most cited author, followed by IRISON (12), 05-68 and GLO-04 (10); JOURNAL OF NEURO-ONCOLOGY (22) was the journal with the largest number of published publications (78), followed by FRONTIERS IN ONCOLOGY (19) and

Frontiers in Oncology | 01 | frontiersin.org

TTFields therapy activates inflammasomes to induce adjuvant immunity in GBM

FOR MORE INFORMATION, USE THE QR CODE:



The Journal of Clinical Investigation RESEARCH ARTICLE

Tumor Treating Fields dually activate STING and AIM2 inflammasomes to induce adjuvant immunity in glioblastoma

Dongling Chen,¹ Song B. Lak,¹ Tarun E. Muthukrishnan,¹ Anil A. Alexandru Calinescu,¹ Matthew Sebastian,¹ Dan Jin,¹ Tianyi Liu,¹ Ashish Choudhury,¹ Mayayeh Rahmani,¹ and David G. Vignani¹

¹Section of Neuro-Oncology and Prostate Cancer, Dana-Farber Cancer Institute, Boston, Massachusetts; and ²Department of Neuro-Oncology, University of Michigan, Ann Arbor, Michigan, USA

Tumor Treating Fields (TTFields), an approved therapy for glioblastoma (GBM) and malignant mesothelioma, employ noninvasive application of low-intensity, intermediate-frequency, alternating electric fields to disrupt the mitotic spindle, leading to chromosome missegregation and apoptosis. Emerging evidence suggests that TTFields may also induce inflammation. However, the mechanism underlying this property and whether it can be harnessed therapeutically are unclear. Here, we report that TTFields induced focal disruption of the nuclear envelope, leading to cytosolic release of long-microsome clusters that intensely recruited and activated 2 major DNA sensors – cyclic-GMP-AMP synthase (cGAS) and AIM2 – and their respective inflammasomes to produce proinflammatory cytokines, type 1 interferons (IFN1s), and IFN1-responsive genes. In lymphoid, murine GBM models, TTFields-treated GBM cells induced antitumor immunity directly and a core set of 42 genes, including a STING- and AIM2-dependent subset. Using single-cell and bulk RNA sequencing of peripheral blood mononuclear cells, we detected robust post-TTFields activation of adaptive immunity in patients with GBM. In PDX-based systems we identified a gene panel signature of TTFields effects on T cell activation and clonal expansion. Collectively, these studies defined a therapeutic strategy using TTFields as cancer immunotherapy in GBM and potentially other solid tumors.

Introduction
Glioblastoma (GBM) is the most common and lethal brain cancer in adults and one of the least immunogenic tumors (1). Recent work has revealed striking immune dysregulation and functional impairment in patients with GBM, including systemic T lymphopenia and anergy and dysfunction of cytokine-producing effector cells (2). GBM tissues also possess a profoundly immunosuppressive or cold tumor microenvironment (TME), characterized by scarce tumor-infiltrating lymphocytes (TILs) and an abundance of inhibitory cells, including myeloid-derived suppressor cells (MDSCs) and regulatory T cells (Tregs). The cold GBM TME expresses high levels of immune checkpoint proteins (3), and is further complicated by tumor cells' production of immunosuppressive factors. In addition, the blood brain barrier (BBB) prevents entry of tumor-associated

antigens to immune cells and also acts, severely hindering immunotherapeutic efforts (4). Overcoming these barriers presents a long-standing, multifaceted, immunotherapy-resistant conundrum. To "heat up" the cold GBM TME, recent efforts have focused on tumor cell-mediated pathways with mixed results, such as dendritic cell-based DC-based vaccines, immune checkpoint blockade, reinvigorating cytokine release, or disrupting BBB integrity to recruit tumor-specific cytotoxic T lymphocytes (CTLs) (5). However, it remains a challenge to engage a direct, core set of tumor cells in reversing the immunosuppressive state of the GBM TME.


In targeting the metabolic, apoptotic, and assembly of macromolecules required for the mitotic spindle structure during mitosis and the associated kinases, dytubulin, tubulin, and cyclin of the cell cycle, Tumor Treating Fields (TTFields) cause chromosome missegregation and aneuploidy and complete cytogenetic separation, respectively, leading to mitotic catastrophe and cell death (6,7). In preclinical systems (8,9), TTFields have also been demonstrated to target the DNA damage repair and breast cancer 1 (BRCA1)-mediated homologous recombination repair pathways by interfering with DNA fork replication (8,10) and induce mitochondrial dysfunction during mitosis (11). Higher order chromosome-associated proteins, kinase-dependent phosphorylation, through structural stabilization of protein-ligand 1 and 2 (L1/L2CA, 9) in breast L1/L2CA (9,11). Recent reports also revealed that TTFields' ability to disrupt the plasma membrane of GBM cells, affecting particles up to 10 nm in size

there is early evidence of efficacy in newly diagnosed GBM patients when TTFields therapy is added to immune checkpoint inhibitors

FOR MORE INFORMATION, USE THE QR CODE:



Overall Survival 2-THE-TOP single arm study vs. external controls



WFNOS 2022 Top 10 Session 2 / March 26 (Sat), 10:15-11:30

BTRT

Abstract
Phase 2 study of pembrolizumab plus TTFields plus temozolomide in patients with newly diagnosed glioblastoma (2-THE-TOP)

David Tran, Ashley Ghaseddin, Dongjiao Chen, Maryam Rahmani
Department of Neurosurgery, Division of Neuro-Oncology, University of Florida, United States

Background: Emerging data indicate that TTFields, the new anti-mitotic treatment for GBM, stimulate immunity via the type-1 interferon (IFN1) pathway of STING and ASK2 inflammasomes. Thus, we hypothesize that TTFields synergize with immune checkpoint inhibitors to induce anti-tumor immunity in GBM.

Methods: We conducted a phase 2 study combining pembrolizumab, TTFields and maintenance TMZ in 26 patients with newly diagnosed GBM (ncGBM). To distinguish immune effects of TTFields from pembrolizumab, TTFields was started at cycle 1 of TMZ and pembrolizumab (200 mg Q2 weeks) at cycle 2. The primary endpoint was PFS vs. the historical control of TTFields plus TMZ (JAMA 318:2306-2316) and immune signatures of TTFields and pembrolizumab by single-cell genomics of PBMCs. Secondary endpoints included toxicity and OS.


Results: As of 09/24/2021, 26 patients with a median age of 60 years were enrolled. Fourteen (54%) had biopsy only or subtotal resection. Nineteen (73%) had immunohistochemical (IHC) and 3 (12%) had an IHC mutation. The median follow-up was 18 and 14.2 months for PFS and OS, respectively. Thirteen (50%) were progression free and 16 (62%) were alive. Of 22 patients with follow-up of months, the median PFS was 3.11 vs. 6.7 months in the control. Six (23%) patients with measurable tumors have achieved partial to complete objective response. We saw opened PD-L1 PD-L2 PD-L3 in 12 patients before and after TTFields and detected robust post-TTFields T cell activation in 11 of 12 patients via the IFN1 regulatory which was strongly correlated with TCRβ clonal expansion (Spearman coefficient $r=0.8$, $p=0.004$). Importantly, we defined a T cell-based gene signature of TTFields effects on TCRβ clonal expansion. The most common adverse events were thrombocytopenia, asthenia, and metabolic disturbances in 4 (15%), 3 (12%), and 2 (8%) patients, respectively.

Conclusions: The triple combination is well tolerated and shows early evidence of efficacy in ncGBM patients. Survival and molecular data will be updated.

Keywords: TTFields immunotherapy; pembrolizumab; STING; single cell analysis

STELLAR Phase 2 trial evaluated TTFields therapy + pemetrexed and cisplatin or carboplatin in MPM

N = 80



Start date: February 2015
Primary completion: April 2018
Study completion: April 2018
Study sites: 13 (Europe)

Primary endpoints:


- OS

Secondary endpoints:

- PFS, ORR (modified RECIST criteria for MPM), safety

MPM patients who used Optune Lua first line achieved 18.2 months median OS

FOR MORE INFORMATION, USE THE QR CODE:



Adapted from Ceresoli GL et al. 2019

Articles

Tumour Treating Fields in combination with pemetrexed and cisplatin or carboplatin as first-line treatment for unresectable malignant pleural mesothelioma (STELLAR): a multicentre, single-arm phase 2 trial

Summary
Background: Tumour Treating Fields (TTFields) are a regional, antineoplastic treatment for solid tumours, which is based on the delivery of low-intensity alternating electric fields. The aim of the STELLAR study was to test the activity of TTFields adjuvant to the therapy in combination with systemic chemotherapy for the first-line treatment of patients with unresectable malignant pleural mesothelioma.

Methods: STELLAR was a prospective, single-arm, phase 2 trial done at 12 European academic and non-academic sites from 16 sites. From 17 February 2016 to 15 October 2018, 80 patients were included in the study. The primary endpoint was overall survival (OS) at 12 months. Secondary endpoints were median OS, 1-year OS rate, and 2-year OS rate. The most common grade 1 or worse adverse events were neutropenia (50%), fatigue (46%), and diarrhoea (30%). The most common grade 2 or worse adverse events were neutropenia (15%), fatigue (10%), and diarrhoea (10%). The most common grade 3 or worse adverse events were neutropenia (10%), fatigue (10%), and diarrhoea (10%). The most common grade 4 or worse adverse events were neutropenia (10%), fatigue (10%), and diarrhoea (10%).

Interpretation: The trial showed encouraging overall survival results, with no increase in systemic toxicity. TTFields did not differ from the therapy combination with pemetrexed and platinum as an active and safe combination for the first-line treatment of unresectable malignant pleural mesothelioma. Further investigations in a randomised trial is warranted.

Funding: Novocure.


Copyright: © 2019 Elsevier Ltd. All rights reserved.

Introduction
Malignant pleural mesothelioma is a rare cancer, characterised by a slow but steady increase in incidence, with an overall survival rate of approximately 12 months. In Europe, there are approximately 4000 new cases diagnosed annually, but it is still increasing in Europe and other countries. Because of the limited progression, most patients are diagnosed with a diffuse disease, inoperable in most cases. The median survival time for patients with unresectable malignant pleural mesothelioma is around 12 months, with platinum and an anti-PD-1 combination being the accepted standard of care since 2017.¹ Carboplatin has had similar efficacy and objective response as cisplatin, but has had a better toxicity profile and more of a maintenance effect with adverse response to administration, and has shown benefit in combination with platinum in the treatment of diffuse malignant mesothelioma in elderly patients.²

The addition of tumour treating fields (TTFields) has been reported to increase survival time by more than 2 months in six reported high-quality randomised, controlled trials.^{3–8} Tumour Treating Fields (TTFields) are a non-invasive, regional, antineoplastic treatment for solid tumours, that is based on the delivery of low-intensity alternating electric fields.

References:
1. Ceresoli GL, et al. *Lancet Oncol*. 2019;20(12):1702–1709.
2. Ceresoli GL, et al. *Lancet Oncol*. 2019;20(12):1702–1709.
3. Ceresoli GL, et al. *Lancet Oncol*. 2019;20(12):1702–1709.
4. Ceresoli GL, et al. *Lancet Oncol*. 2019;20(12):1702–1709.
5. Ceresoli GL, et al. *Lancet Oncol*. 2019;20(12):1702–1709.
6. Ceresoli GL, et al. *Lancet Oncol*. 2019;20(12):1702–1709.
7. Ceresoli GL, et al. *Lancet Oncol*. 2019;20(12):1702–1709.
8. Ceresoli GL, et al. *Lancet Oncol*. 2019;20(12):1702–1709.

EF-15 Phase 2 trial evaluated TTFields therapy + pemetrexed in NSCLC



Start date: May 2008

Primary completion: July 2011

Study completion: July 2011

Study sites: 4 (Switzerland)

Primary endpoints:


- Device related toxicity (P1), Time to in-field progression (P2)

Secondary endpoints:

- OS, ORR, time to systemic progression, safety

TTFields therapy together with pemetrexed improved disease control within the treatment field in second line NSCLC

FOR MORE INFORMATION, USE THE QR CODE:



Median in-field PFS Median PFS Median OS 1yr Survival

TTFields + Pemetrexed

6.5 mo

5.0 mo

13.8 mo

57.0%

Pemetrexed alone

n/a


2.9 mo

8.3 mo

29.7%



LUNAR Phase 3 trial evaluated TTFields therapy + SOC in metastatic NSCLC, post-platinum



Start date: December 2016

Primary completion: December 2022

Study completion: December 2022

Study sites: 124

Primary endpoints:

- OS

Secondary endpoints:


- OS (by cohort), PFS, ORR, QoL, safety

TTFields therapy together with either standard of care therapies or immune checkpoint inhibitor improved overall survival in second-line NSCLC



FOR MORE INFORMATION, USE THE QR CODE:




Overall survival (ITT population)



Overall survival (ICI-treated patients)



METIS Phase 3 trial evaluated TTFields therapy + supportive care in NSCLC brain metastases, following SRS



Start date: October 2016
Primary completion: March 2023
Study sites: 125

Primary endpoints:

- Time to intracranial progression

Secondary endpoints:


- Time to neurocognitive failure, OS, radiological response rate, time to 2nd intracranial progression, time to 1st and 2nd progression by cohort (1-4 metastases, 5-10 metastases), rate of intracranial progression at two-month intervals, time to distant progression, rate of cognitive decline, neurocognitive failure-free survival, quality of life, adverse events

TTFields therapy with supportive care following SRS improved time to intracranial progression in patients with brain metastases from NSCLC

FOR MORE INFORMATION, USE THE QR CODE:



Time to Intracranial Progression



CENTRAL NERVOUS SYSTEM TUMORS

2008

Check for updates


David Albertson, David

Results from METIS (EF-25), an international, multicenter phase III randomized study evaluating the efficacy and safety of tumor treating fields (TTFields) therapy in NSCLC patients with brain metastases

Mehta P, Mehta N, Smith-Horowitz S, et al. *J Clin Oncol* 2024;42:16 suppl.2008

Background: For non-small cell lung cancer (NSCLC) with brain metastases (BM), stereotactic radiosurgery (SRS) is the current preferred therapy. Due to frequent intracranial failures, there is a high unmet need for salvage therapies. Whole brain radiotherapy (WBRT) reduces intracranial failure but is often limited by cognitive consequences. Tumor Treating Fields (TTFields), an electric field-based therapy, has been shown to improve survival and safety in patients with glioblastoma and metastatic NSCLC. Phase 3 METIS trial (NCT03759148) was designed to evaluate the efficacy and safety of TTFields therapy in NSCLC patients with BM treated with SRS, specifically on terms of lengthening time to intracranial progression without cognitive decline. Methods: Metastatic negative (M-) NSCLC patients with 1-10 BM were randomized 1:1 to receive stereotactic radiosurgery (SRS) followed by Tumor Treating Fields (TTFields), SRS followed by best supportive care (BSC) or SRS followed by BSC. Patients with Karnofsky Performance Status (KPS) \geq 70, newly diagnosed with no resectable or \geq 10 mm supratentorial brain metastases suitable for SRS and receiving optimal intracranial disease therapy were included. Exclusions were prior WBRT and single operable or recurrent brain metastases. Primary endpoint was time to first intracranial progression (M0+M1) based on cumulative risk. Patients were followed every two months until second intracranial progression. Cognitive and patient quality of life (QoL) were evaluated. Results: Between July 2017 and September 2022, 289 patients were randomized. Baseline characteristics were balanced: median age was 63.3 (range 17-84) years, 33.6% female, majority of patients had a KPS \geq 80, median time from initial NSCLC diagnosis was 1.8 months (range 0.2-15.1), 77% had adenocarcinoma. Median treatment duration of TTFields was 16 weeks, with median usage time of 67%. Primary endpoint, time to intracranial progression from SRS, was significantly prolonged with SRS followed by TTFields therapy with BSC vs. TTFields plus BSC arm (median of 21.9 vs. 11.3 months); HR=0.67 (95% CI 0.46-1.0), $p=0.016$. TTFields-related AEs were mainly dermatologic, and Grade \geq 3 TTFields therapy also improved deterioration-free survival of global health status, physical functioning, and fatigue, according to QoL, and did not negatively impact cognition. Conclusions: METIS study met its primary endpoint, demonstrating that TTFields therapy following SRS significantly improved NSCLC patients with BM, significantly prolonging time to intracranial progression and could postpone WBRT, without QoL and cognitive decline. Clinical trial information: NCT03759148. Research sponsor: Novocure GmbH.

PANOVA phase 2 trial evaluated TTFields therapy + gemcitabine +/- nab-paclitaxel in pancreatic cancer



Start date: Nov 2013

Primary completion date: Dec 2017

Study completion date: Dec 2017

Study sites: 6 (Europe)

Primary endpoint:


- Safety

Secondary endpoints:

- TTFields monthly usage, PFS, OS


TTFields therapy together with chemotherapy were well tolerated for patients with advanced pancreatic cancer

FOR MORE INFORMATION, USE THE QR CODE:




OS	
Median, mo	NR
95% CI	8.4–NA
1-year survival	72%

	Median PFS	Median OS	One-year Survival	Partial Response Rate	Stable Disease
TTFields + gemcitabine	8.3 mo	14.9 mo	55%	30%	30%
gemcitabine alone	3.7 mo	6.7 mo	22%	7%	28%
TTFields + gemcitabine + nab-paclitaxel	12.7 mo	Not yet reached	72%	40%	47%
gemcitabine + nab-paclitaxel alone	5.5 mo	8.5 mo	35%	23%	27%



PANOVA-3 Phase 3 trial evaluated TTFields therapy + gemcitabine + nab-paclitaxel in unresectable, locally advanced pancreatic cancer



Start date: February 2018
Primary completion: October 2024
Study completion: October 2024
Study sites: 199

Primary endpoints:

- OS

Secondary endpoints:

- PFS, local PFS, ORR, 1-year survival rate, QoL, pain-free survival, puncture-free survival, resectability rate, safety


MET PRIMARY ENDPOINT

encouraging response rate and durability signals in EF-31 phase 2 gastric cancer trial

FOR MORE
INFORMATION, USE
THE QR CODE:



EF-31 PHASE 2 PILOT TRIAL DESIGN¹




	OBJECTIVE RESPONSE RATE	MEDIAN PROGRESSION-FREE SURVIVAL	DURATION OF RESPONSE	ONE-YEAR SURVIVAL
TTFields + chemotherapy	50%	7.8mo	10.3mo	72%
SOC chemotherapy ²	41-45%	6.9mo	6.9mo	48%

encouraging signals in liver cancer despite poor prognosis and low treatment exposure in HEPANOVA phase 2 trial

FOR MORE INFORMATION, USE THE QR CODE:



HEPANOVA PHASE 2 PILOT TRIAL DESIGN²



76%

DISEASE CONTROL RATE
(n=21)

VS. 43% CONTROL³

9.5%

OBJECTIVE RESPONSE RATE
(n=21)

VS. 4.5% CONTROL

91%

DISEASE CONTROL RATE

18%

OBJECTIVE RESPONSE RATE

patients that received ≥ 12 wks of TTFields (n=11)

2025-2026 anticipated clinical development milestones


	TRIAL	TTFIELDS +	PHASE 2	PHASE 3	APPROVED
CNS indications	glioblastoma	EF-14	TMZ	APPROVED ✓	
		TRIDENT	TMZ + radiation	DATA IN 2026	
		KEYNOTE D58	TMZ + pembrolizumab	enrolling	
	brain metastases	EF-11	monotherapy (recurrent GBM)	APPROVED ✓	
		METIS	monotherapy	SUBMISSION IN 2025	
torso indications	non-small cell lung cancer	LUNAR	docetaxel or PD-L1 inhibitor (2L)	APPROVED ✓	
		LUNAR-2	pembrolizumab + platinum (1L)	enrolling	
		LUNAR-4	pembrolizumab (2L retreatment)	enrolling	
	mesothelioma	STELLAR	pemetrexed + cisplatin/carboplatin	APPROVED ✓	
	pancreatic cancer	PANOVA-3	nab-paclitaxel + gemcitabine (LAPC)	SUBMISSION IN 2025	
		PANOVA-4	atezolizumab + nab-paclitaxel + gemcitabine (MPC)	DATA IN 2026	



real-world and
routine clinical care
study evidence

real-world evidence showed ndGBM median overall survival extension by over 12 months in the high use TTFields group

FOR MORE INFORMATION, USE THE QR CODE:



Neuro-Oncology Advances

4715-11-8-2022 | <https://doi.org/10.1093/advances/nwac012> | Advance Access date 19 September 2022

Determinants of tumor treating field usage in patients with primary glioblastoma: A single institutional experience

Matthew T. Ballo, Kathleen W. Qualls, L. Madison Michael, Jeffrey M. Sonenson, Brandon Baughman, Saradeen Karri-Walsh, and Manjari Pandey

Department of Radiation Oncology, West Cancer Center & Research Institute, Memphis, Tennessee, USA (M.T.B., K.W.Q.); Neurosurgery, Sitemp Memphis Neurosurgery Clinic, Memphis, Tennessee, USA (L.M.M., J.M.S., B.B.); Department of Medical Oncology, West Cancer Center and Research Institute, Memphis, Tennessee, USA (S.K.W., M.P.)

Corresponding Author: Matthew T. Ballo, MD, Department of Radiation Oncology, West Cancer Center and Research Institute, 7985 Wolf River Blvd, Germantown, TN 38138, USA (mailto:mballo@wcri.com)

Abstract
Background. Determinants of tumor treating fields (TTFields) usage in patients receiving combined modality therapy for primary CNS wild-type glioblastoma are currently unclear.
Methods. Ninety-one patients underwent maximal debulking surgical resection, completed external beam radiotherapy with concurrent temozolomide (TMZ), and initiated adjuvant TMZ with or without TTFields. We performed a retrospective analysis of patient, tumor, and treatment-related factors that affected TTFields usage.
Results. We identified three TTFields usage subgroups: 32 patients that declined TTFields, 40 patients that started, but had monthly compliance of less than 75% or used it for less than 2 months, and 19 patients who used TTFields for 2 or more months and maintained average monthly compliance greater than 75%. With 26.5 months median follow-up for surviving patients, the 1- and 3-year actuarial overall survival for all patients was 80% and 58%, respectively. On multivariate analysis TTFields use (P = .03), extent of surgical resection (P = 0.02), and MGMT methylation status (P = .01) were significantly associated with overall survival. TTFields usage was explored as a continuous variable and higher average usage was associated with longer overall survival (P = .03). There was no relationship between patient, tumor, or treatment-related factors and a patient's decision to use TTFields.
Conclusions. No subgroup of patients was more or less likely to initiate TTFields therapy and no subgroup was more or less likely to use TTFields as prescribed. The degree of TTFields compliance may be associated with improved survival independent of other factors.


Key Point
 • It is reasonable to offer all patients with primary glioblastoma TTFields therapy as we could not identify a group that was more or less likely to discontinue therapy or unable to initiate therapy. Patients benefit from TTFields regardless of tumor or patient characteristics.

Glioblastoma is the most common and aggressive primary malignant brain tumor diagnosed in adults and has a poor survival has historically been only 14.6 months.^{1,2} Tumor treating fields (TTFields) represent a novel therapy in the treatment of glioblastoma. TTFields deliver low-intensity, intermediate-frequency (200 kHz) alternating electric fields

© The Author(s) 2022. Published by Oxford University Press, the Society for Neuro-Oncology and the European Association of Neuro-Oncology. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted reuse, distribution, and reproduction in any medium, provided the original work is properly cited.

real-world evidence validates EF-14 with statistically significant improvement in PFS and OS in Chinese patients with ndGBM

FOR MORE INFORMATION, USE THE QR CODE:



	Median OS	Median PFS
TTFields + TMZ (n=63)	21.8 mo	16.0 mo
TMZ alone (n=204)	15.0 mo	11.0 mo

Journal of Clinical Medicine MDPI

Article
Tumor Treating Fields Combine with Temozolomide for Newly Diagnosed Glioblastoma: A Retrospective Analysis of Chinese Patients in a Single Center

Chenxi Chen ^{1,2,3,4,5,6,7,8,9,10,11}, Han Yu ^{1,2,3,4,5,6,7,8,9,10,11}, Kun Song ^{1,2,3,4,5,6,7,8,9,10,11}, Yi Zhang ^{1,2,3,4,5,6,7,8,9,10,11}, Jianyan Zhang ^{6,7,8,9,10,11}, Yang Wang ⁴, Xiaoting Sheng ⁴, Lingchen Chen ^{1,2,3,4,5,6,7,8,9,10,11} and Zhiyong Qiu ^{1,2,3,4,5,6,7,8,9,10,11}

¹ Department of Neurosurgery, Shanghai Hospital, Fudan University, Shanghai 20000, China; ² National Center for Neurological Disorders, Shanghai 20000, China; ³ Shanghai Key Laboratory of Brain Function Restoration and Neural Regeneration, Shanghai 20000, China; ⁴ Neurological Institute of Fudan University, Shanghai 20000, China; ⁵ Shanghai Clinical Medical Center of Neurosurgery, Shanghai 20000, China; ⁶ Institute of Clinical Neurophysiology, Shanghai JITPC, China; ⁷ Branch of Clinical Epidemiology and Evidence-Based Medicine, Shanghai Medical Association, Shanghai 20000, China; ⁸ Department of Radiation Oncology, Shanghai Hospital, Fudan University, Shanghai 20000, China; ⁹ Correspondence: chenxi@shsmu.edu.cn (C.C.) or zhiyongqiu@shsmu.edu.cn (Z.Q.); ¹⁰ These authors contributed equally to this work; ¹¹ These authors contributed equally to this work

Abstract: Introduction: TTFields plus temozolomide (TTFields/TMZ) extended survival versus TMZ alone in newly diagnosed glioblastoma (GBM) patients in the EF-14 trial. We have reported a retrospective analysis of newly diagnosed Chinese GBM patients who received TTFields/TMZ treatment and TMZ treatment from August 2019 to May 2022 in Shanghai Hospital in Shanghai, Mainland, Overall survival (OS) and progression-free survival (PFS) curves were constructed using the Kaplan-Meier method. A Cox proportional hazards regression model, propensity score matched data, and inverse probability of treatment weighting (IPTW) based on propensity score were used to assess the effect of TTFields and account for confounding factors. Results: In the preliminary analysis, the median PFS in TTFields/TMZ group was 16 months (95% CI: 9.4–24.0) versus 11 months (95% CI: 6.52 to 15.62) in TMZ group ($p < 0.05$). Median overall survival was 21.8 months (95% CI: 17.4–NA) with TTFields/TMZ versus 15 months (95% CI: 11.48–19% CI: 15–30) with TMZ alone. The multivariate analysis identified younger age, SEEP scheme, EBH status, and TTFields as an overall prognostic factors. After PSM adjustment, the results among the groups were similar, except that the multivariate use of TMZ/TMZ predictor remained high in the TMZ group (12/15 months; $p < 0.001$). Upon IPTW survival analysis, TTFields was associated with a significantly lower risk of death (HR = 0.39 in OS, 95% CI: 0.19–0.81) and progression (HR = 0.20, 95% CI: 0.11–0.40) compared with TMZ group. Conclusion: In the final analysis of our single-center Chinese patients with glioblastoma, adding TTFields to temozolomide chemotherapy resulted in statistically significant improvement in PFS and OS.


Keywords: glioblastoma; tumor treating fields; chemotherapy; retrospective cohort

Copyright: © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<http://creativecommons.org/licenses/by/4.0/>).

J. Clin. Med. **2023**, *11*, 5855. <https://doi.org/10.3390/jcm11195855> <https://www.mdpi.com/journal/jcm>

post-approval study supports safety and efficacy profile of TTFields therapy in ndGBM Japanese patients, validating EF-14 improved survival rates

FOR MORE INFORMATION, USE THE QR CODE:



	1-year survival	2-year survival
TMZ alone ¹	65%	31%
TTFields + TMZ (n=14)	77.9%	53.6%
TTFields + TMZ high use group (n=21)	89.5%	65.8%

JJCO Japanese Journal of Clinical Oncology

Volume 52, Number 14, December 2023
 ISSN 1340-3445
<https://doi.org/10.1093/jjco/tyad001>

Original Article

Safety and efficacy of tumour-treating fields (TTFields) therapy for newly diagnosed glioblastoma in Japanese patients using the Novo-TTF System: a prospective post-approval study

Ryo Nishikawa^{1,2}, Fumiyuki Yamasaki³, Yoshiki Arakawa⁴, Yoshihiro Muragaki⁵, Yoshitaka Naita⁶, Shota Tanaka⁷, Shigeru Yamaguchi⁸, Akitsue Mukasa⁹ and Masayuki Kanomura¹⁰

¹Department of Neuro-Oncology/Neurosurgery, Satama Medical University International Medical Center, Satama, Japan; ²Department of Neurosurgery, Hiroshima University Hospital, Hiroshima, Japan; ³Department of Neurosurgery, Kyoto University Graduate School of Medicine, Kyoto, Japan; ⁴Department of Neurosurgery, Tokyo Women's Medical University Hospital, Tokyo, Japan; ⁵Department of Neurosurgery and Neuro-Oncology, National Cancer Center Hospital, Tokyo, Japan; ⁶Department of Neurosurgery, The University of Tokyo Hospital, Tokyo, Japan; ⁷Department of Neurosurgery, Hokkaido University Hospital, Sapporo, Japan; ⁸Department of Neurosurgery, Kumamoto University Hospital, Kumamoto, Japan and ⁹Department of Neurosurgery, Shinshu University Hospital, Sendai, Japan


¹⁰For reprints and all correspondence: Ryo Nishikawa, Satama Medical University International Medical Center, Department of Neuro-Oncology/Neurosurgery, 1300-1 Satama, Hokkaido-shi, Satama 069-1388, Japan. E-mail: rnoh31205@gmail.com

Received 16 July 2023; Revised 17 November 2023; Editorial Decision 28 December 2023; Accepted 2 January 2024


Abstract
 Background: Tumour-treating fields therapy is a locoregional, anti-cancer treatment. Efficacy and safety of tumour-treating fields therapy in adults with newly diagnosed glioblastoma were demonstrated in the pivotal phase 3 EF-14 study (NCT02102010). Here, we report post-approval data of tumour-treating fields therapy in Japanese patients with newly diagnosed glioblastoma. Methods: Unpublished post-marketing surveillance data from Japanese patients with newly diagnosed glioblastoma treated with tumour-treating fields therapy (December 2016–June 2020) were retrospectively analysed. The primary endpoints were skin, neurological and psychiatric adverse events. The secondary endpoints were 1- and 2-year overall survival rates, and the 6-month progression-free survival. Adverse events were analysed using MedDRA v24.0. The overall survival and progression-free survival were assessed using the Kaplan–Meier survival analysis (log-rank testing). The Cox proportional hazard regression analyses were also performed. Results: Forty patients with newly diagnosed glioblastoma were enrolled (85.2% male; median age 59 years; median baseline Karnofsky Performance Scale score 90). The most common tumour-treating fields-therapy-related adverse event was beneath-army local skin reaction (100% of patients). The adverse events were mostly mild to moderate in severity. Neurological disorders were observed in 2.5% patients (one patient reported dysesthesia). No psychiatric disorders were reported. The 1- and 2-year overall survival rates were 77.9% (95% CI 65.8–83.2) and 53.6% (95% CI 46.7–62.7), respectively. The 6-month progression-free survival was 77.5% (95% CI 67.8–83.2). These survival rates were similar to those reported in the EF-14 study. Conclusion: TTFields therapy was well tolerated and showed promising efficacy in Japanese patients with newly diagnosed glioblastoma. These results support the use of TTFields therapy in Japanese patients with newly diagnosed glioblastoma.

long term study of ndGBM patients, covering 18-year period, confirms TTFields' positive effect on PFS and OS

FOR MORE INFORMATION, USE THE QR CODE:



	median PFS	median OS
TTFields + TMZ (n=55)	19.75 mo	31.67 mo
TMZ alone ¹	12.45 mo	24.80 mo




1) Results from the EF-14 study, Stupp R et al. JAMA. 2017;318(23):2306–2316. ndGBM, newly diagnosed glioblastoma; OS, median overall survival; PFS, median progression-free survival; TMZ, temozolomide; TTFields, Tumor Treating Fields. Mild to moderate skin irritation was the most common device related side effect. Vymazal J, et al. Front. Oncol. 12:1014455. doi: 10.3389/fonc.2022.1014455

TIGER study of routine clinical care in German ndGBM patients corroborates overall survival and safety outcomes from EF-14


FOR MORE INFORMATION, USE THE QR CODE:



Overall Survival



Progression-Free Survival



Tumor Treating Fields therapy in patients with newly diagnosed glioblastoma: Long-term survival results from TTFields in Germany in routine clinical care (TIGER) study
 Oliver Bähr¹, Chahab Tahakolat², Bastian Falkau³, Roland Goldammer⁴, Martin Glat⁵

Background

- TTFields therapy is a non-invasive, non-drug treatment for glioblastoma (GBM) that has been shown to improve overall survival (OS) in patients with newly diagnosed GBM (ndGBM) in the EF-14 phase III trial.
- The TIGER study is a long-term, retrospective analysis of patients with ndGBM who received TTFields therapy in routine clinical care in Germany.
- The primary objective of the TIGER study is to evaluate the long-term survival outcomes of patients with ndGBM who received TTFields therapy in routine clinical care.
- The secondary objectives of the TIGER study are to evaluate the safety, quality of life, and health economics of TTFields therapy in routine clinical care.

Methods

The TIGER study included patients with ndGBM who received TTFields therapy in routine clinical care in Germany between 2010 and 2020. The study included patients who received TTFields therapy as a first-line treatment for ndGBM.

Results

The median OS of patients with ndGBM who received TTFields therapy in routine clinical care was 19.6 months (95% CI: 17.9 - 22.4). The median PFS was 10.2 months (95% CI: 9.4 - 11.4).

Conclusions

- TIGER is the largest prospective study to date on routine clinical practice in GBM.
- TIGER confirms TTFields therapy as a safe, effective, and well-tolerated treatment for ndGBM.
- Patients with TTFields therapy demonstrate promising long-term survival outcomes in patients with ndGBM.

review article identifies TTFields therapy as one of few factors driving increased overall survival in GBM patients since the 2005 Stupp-protocol

FOR MORE INFORMATION, USE THE QR CODE:




Neff et al.

Characteristic	N	HR ¹	95% CI ¹	p-value
Age (years)	19,414	1.02	1.02, 1.03	<0.001
Sex				
Female	8,046	—	—	reference
Male	11,368	1.10	1.07, 1.14	<0.001
Elkhauser Comorbidity Score	19,414	1.01	1.01, 1.01	<0.001
Tumor-Treating Fields (ever)				
No	16,353	—	—	reference
Yes	3,061	0.77	0.73, 0.80	<0.001
Received radiation or radiosurgery (ever)				
No	7,370	—	—	reference
Yes	12,044	0.88	0.85, 0.91	<0.001
Bevacizumab (ever)				
No	15,741	—	—	reference
Yes	3,673	0.85	0.82, 0.88	<0.001

In this commercially insured dataset, TTFields improved OS to a greater extent (HR=0.77) vs. Bevacizumab (HR=0.85) or Radiation use (HR=0.88)

TTFields subset n=3,061 over 6 years

Mrugala et al.



AEs were consistent with the safety profile from the pivotal EF-11 and EF-14 clinical studies

n=23,822 over 11 years





tumor treating fields mechanism of action appendix

patients with aggressive solid tumors often face suboptimal survival outcomes, despite advancements in treatment modalities

These outcomes are due to diverse treatment challenges, including:



Therapeutic tumor resistance



Drug-to-drug interactions




Additive systemic toxicities

With a poor survival outlook, physicians and patients need additional treatment strategies

Tumor Treating Fields (TTFields) are electric fields that exert physical forces to kill cancer cells via a variety of mechanisms



TTFields spare healthy cells because they have different properties than cancer cells across a range of tumor types



a growing body of evidence supporting multiple mechanisms of action

FOR MORE
INFORMATION, USE
THE QR CODE:



- Preclinical research has shown interference with cancer cell motility and migration, activation of anti-tumor immunity, downregulation of genes important for DNA damage repair, and other potential mechanisms
- May demonstrate enhanced effects across solid tumor types when used with chemotherapy, radiotherapy, immune checkpoint inhibition, or PARP inhibition in preclinical models

JOURNAL ARTICLE ACCEPTED MANUSCRIPT

Anti-cancer mechanisms of action of therapeutic alternating electric fields (tumor treating fields [TTFields])

Shadi Shams, Chirag B Patel


Journal of Molecular Cell Biology, mjac047, <https://doi.org/10.1093/jmcb/mjac047>
Published: 16 August 2022

PDF Split View Cite Permissions Share

Abstract



Despite improved survival outcomes across many cancer types, the prognosis remains grim for certain solid organ cancers including glioblastoma and pancreatic cancer. Invariably in these cancers, the control achieved by time-limited interventions such as traditional surgical resection, radiation therapy, and chemotherapy is short-lived. A new form of anti-cancer therapy called therapeutic alternating electric fields (AEFs) or tumor treating fields (TTFields) has been shown, either by itself or in combination with chemotherapy, to have anti-cancer effects that translate to improved survival outcomes in patients. Although the pre-clinical and clinical data are promising, the mechanisms of TTFields are not fully elucidated. Many investigations are underway to better understand how and why TTFields is able to selectively kill cancer cells and impede their proliferation. The purpose of this review is to summarize and discuss the reported mechanisms of action of TTFields from pre-clinical studies (both *in vitro* and *in vivo*). An improved understanding of how TTFields works will guide strategies focused on the timing and combination of TTFields with other therapies, to further improve survival outcomes in patients with solid organ cancers.

Tumor Treating Fields have multiple, distinct mechanisms of action




TTFields have been shown to disrupt mitosis in cancer cells by exerting physical forces on their polar components

FOR MORE INFORMATION, USE THE QR CODE:



TTFields have been shown to alter the organization and dynamics of the cytoskeleton, disrupting cancer cell motility and migration

FOR MORE INFORMATION, USE THE QR CODE:



a) microtubule; b) TTFields; c) tubulin aligned with field; j) actin fiber; k) integrin; l) focal adhesion; m) extracellular matrix.

A model illustrating the mechanism by which TTFields modulates cancer cell motility.

- (1) Microtubules are required to specify the direction of cell movement. GEF-H1 catalytic activity is downregulated through microtubule binding.
- (2) TTFields exert directional forces on polar tubulins leading to their alignment in the direction of the field. This, in turn, leads to the reorganization of the microtubule network resulting in changes in the abundance of microtubules and initiation of the GEF-H1/RhoA/ROCK signaling pathway
- (3) to increase actin bundling
- (4) and formation of focal adhesions,
- (5) which disrupt cell polarity and migration directionality.

cancers

Tumor Treating Fields (TTFields) Hinder Cancer Cell Motility through Regulation of Microtubule and Actin Dynamics

Tali Voloshin¹, Roni Sara Schneiderman¹, Alexandra Volodin, Reuben Rudy Shmitz, Noa Keren, Elmer Zevi, Ilitak Keren, Anat Klein-Goldberg, Ron Tam, Moshe Giladi¹, Zeev Ben-Zur, Uri Weisberg and Nissim Kashi

Novocure Ltd., Super Building, MASHIM Center Haifa 31091, Israel; voloshin@novocure.com (T.V.); schneiderman@novocure.com (R.S.S.); volodin@novocure.com (A.); shmitz@novocure.com (R.S.); keren@novocure.com (E.Z.); ilitak@novocure.com (I.K.); koren@novocure.com (A.K.); zevi@novocure.com (E.Z.); klein-goldberg@novocure.com (A.G.); tam@novocure.com (R.T.); moshe.giladi@novocure.com (M.G.); zeev.ben-zur@novocure.com (Z.B.); uris@novocure.com (U.W.); nissim.kashi@novocure.com (N.K.)

¹ These authors have contributed equally to this work.

Received: 21 September 2020; Accepted: 14 October 2020; Published: 17 October 2020

Simple Summary: Tumor Treating Fields (TTFields), encompassing alternating electric fields within the intermediate frequency range, is an anticancer treatment delivered to the tumor region through transducer arrays placed non-invasively on the skin. Although established as an anti-neoplastic treatment modality, the anti-metastatic potential of TTFields and their effect on signal transduction dynamics during cellular motility warrant further investigation. In this study, we report that TTFields application induces changes in microtubule organization leading to interference with the directionality and abundance of cancer cell migration. We show that these changes in microtubule organization result in activation of GEF-H1/RhoA/ROCK signaling pathway and the consequent formation of focal adhesions and changes in actin cytoskeleton architecture. Together, these results propose a novel mechanism by which TTFields induce changes in microtubule and actin organization and dynamics, thereby disrupting processes important for polarity generation and motility in cancer cells.

Keywords: Tumor Treating Fields (TTFields); microtubule; alternating electric fields; intermediate frequency range; anti-neoplastic treatment


TTFields are homogeneously delivered to the tumor region through 2 pairs of transducer arrays placed on the skin. This novel treatment modality has been FDA-approved for use in patients with glioblastoma and multiple drug resistant metastatic breast cancer. However, the effectiveness, efficacy and safety, and its currently under investigation in other types of solid tumors. TTFields were shown to induce an anti-metastatic effect by exerting bidirectional forces on highly polar intracellular cytoskeleton, such as tubulin and actin microtubules, leading to abnormal microtubule polymerization during spindle formation as well as aberrant cleavage furrow formation. Previous studies have demonstrated that TTFields inhibit metastatic progression in cancer cells. However, the consequences of TTFields application on cytoskeleton dynamics remain undetermined. In this study, methods defined in contribution to study the effects of TTFields on cancer cell motility through regulation of microtubule and actin dynamics included confocal microscopy, computational tools, and biochemical analyses. Mechanisms by which TTFields treatment disrupted cellular polarity were (1) interference with microtubule assembly and disassembly; (2) altered regulation of G-actin microtubule exchange factor (GEF-H1), RhoA family member (RhoA), and Rho-associated coiled-coil kinase (ROCK) activity; and (3) induced formation of radial protrusions of peripheral actin filaments and focal adhesions. Overall, these data identified discrete effects of TTFields that disrupt processes crucial for cancer cell motility.

www.mdpi.com/journal/cancers

Copyright: © 2020 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<http://creativecommons.org/licenses/by/4.0/>).

TTFields-mediated cell disruption activates the immune system and triggers a downstream antitumor cell response

FOR MORE INFORMATION, USE THE QR CODE:



TTFields induces downstream immunogenic cell death, including release of DAMPs (damage-associated molecular patterns)




TTFields downregulate genes important for DNA damage repair

FOR MORE INFORMATION, USE THE QR CODE:



- TTFields disrupt DNA damage repair in cancer cells by downregulating genes that are part of the well-known FA-BRCA pathway^{1,2}



OPEN

Tumor-treating fields elicit a conditional vulnerability to ionizing radiation via the downregulation of BRCA1 signaling and reduced DNA double-strand break repair capacity in non-small cell lung cancer cell lines

Narasimha Kumar Karanam¹, Kalyanram Srinivasan¹, Lianhua Ding¹, Brock Blain¹, Debabrata Saha¹ and Michael D Story^{1,2}

Cell Death and Disease 2017; 8: e2171. doi:10.1038/cddis.2017.136

The use of tumor-treating fields (TTFields) has revolutionized the treatment of recurrent and newly diagnosed glioblastoma (GBM). TTFields are low-frequency, intermediate frequency alternating electric fields that are applied to tumor regions and cells using non-invasive arrays. The predominant mechanism by which TTFields are thought to kill tumor cells is the disruption of mitosis. Using the non-small cell lung cancer (NSCLC) cell lines we found that there is a variable response to cell proliferation and cell killing between the NSCLC cell lines that are independent of cell status. TTFields are most potent in the G2M population, with a concomitant reduction in S-phase cells followed by the appearance of a sub-G1 population indicative of apoptosis. Temporal changes in gene expression during TTFields exposure were evaluated to identify molecular signaling changes underlying the differential TTFields response. The most differentially expressed genes were associated with the cell cycle and cell proliferation pathways. However, the expression of genes found within the BRCA1-DNA damage response were significantly downregulated (P < 0.05) during TTFields treatment. DNA double-strand break (DSB) repair foci increased when cells were exposed to TTFields as did the appearance of ataxia telangiectasia, suggesting an ataxia telangiectasia mechanism responsible for cell death involving DNA repair. Exposing cells to TTFields immediately following ionizing radiation resulted in increased chromosomal aberrations and a reduced capacity to repair DNA DSBs, which were fully reversible in at least a portion of the enhanced cell killing seen with the combination. These findings suggest that TTFields induce a state of BRCA1 loss leading to a conditional susceptibility resulting in enhanced sensitivity to ionizing radiation and provide a strong rationale for the use of TTFields as a combined modality therapy with radiation or other DNA-damaging agents.

Cell Death and Disease (2017) 8, e2171. doi:10.1038/cddis.2017.136. published online 26 March 2017

Lung cancer is the second most prevalent cancer and the leading cause of cancer-related death in the United States.^{1,2} Non-small cell lung cancer (NSCLC) is the most prevalent type, accounting for ~85% of new cases.^{3,4} A plethora of treatment options, each including surgical resection, chemotherapy, radiation therapy and immunotherapy,^{5,6} have emerged over time for patients with stage I-III NSCLC, but the overall 5-year survival remains dismal for most of the options. Systemic cancer care for patients with late-stage NSCLC and for the 15%, 15%, 15% of patients with recurrent, metastatic and refractory NSCLC, respectively, remains limited, highlighting the need for novel treatment modalities that can be utilized alone or in combination with conventional therapies to increase survival rates.

The advent of Tumor-Treating Fields (TTFields), a novel physical treatment modality, has been effective for the treatment of adult, recurrent glioblastoma (GBM) and recurrent tumors.^{7,8} TTFields are non-invasive and deliver a non-thermal, non-ionizing, intermediate frequency (150-200 kHz) alternating electric field across the tumor site.^{9,10} TTFields create a heterogeneous intracellular environment that induces a dielectrophoretic movement of polar molecules towards the region of higher field intensity, effectively preventing cytoskeletal and other critical biochemical functions.¹¹ As such, TTFields prevent adult lung cancer cells through the inhibition of cell proliferation, effectively sparing non-dividing normal cells in addition. TTFields do not ablate nerves and muscle because of their high frequency, and do not penetrate the brain due to their non-invasive nature. The safety and efficacy of TTFields in combination with immunotherapy, radiation therapy, and chemotherapy for the treatment of recurrent and refractory NSCLC is being investigated in clinical trials.¹² Clinical trials are ongoing or ready for initiation in patients with recurrent NSCLC, lung primary and metastatic cancers (www.novocure.com).

TTFields are known to disrupt cellular proliferation and induce apoptotic pathways in dividing cancer cells across a variety of human solid tumor types cell lines.¹³ Proliferation of larger fractions of the mitotic spindle apparatus and the inhibition of other mitotic structures have been proposed as the mechanism by which TTFields kill dividing cells.^{14,15} Specifically, TTFields require viable or non-dividing replication and the mislocalization of spindles. This results in

TTFields is a highly versatile first-in-class treatment modality

FOR MORE
INFORMATION, USE
THE QR CODE:




- TTFields therapy has significant potential for broad applicability across solid tumor types and lines of therapy
- Investigation of TTFields therapy is ongoing across clinical trials in multiple tumor types
- In approved indications, TTFields therapy is well tolerated, suggesting a low risk of additive systemic toxicity when used with other cancer treatment modalities




TTFields therapy can be added to cancer treatment modalities in approved indications


TTFields **demonstrate enhanced effects** across multiple solid tumor types, when used concomitantly with each of the following:




Chemotherapy



Radiation therapy (RT)



Targeted therapies



Immuno-oncologic (IO) agents