



Checklist for starting treatment with Optune Lua®

Optune Lua is a wearable, portable, FDA-approved device used together with PD-1/PD-L1 inhibitors (immunotherapy) or docetaxel. It is indicated for adult patients with metastatic non-small cell lung cancer (mNSCLC) who have progressed on or after a platinum-based regimen.

Use this guide to help you get ready for your Optune Lua treatment.

Actor Portrayal

Once your doctor has prescribed Optune Lua and you are ready to start therapy:

- A **MyNovocure®** team member will contact you for a welcome call
- A **Device Support Specialist (DSS)** will also contact you to schedule an appointment to start therapy. It is recommended to have a family member or friend join you for this appointment if possible. They can learn how to help you with placing and changing your arrays during your therapy
- The **Optune Lua device** can be shipped to your home before your appointment or brought to your home when you start treatment. Your DSS will discuss these options on the start call
- An **electric shaver** can be provided prior to your start appointment, if needed

Make sure you have the following items at home before starting therapy:

- ☐ A **skin barrier film** recommended by your doctor
- ☐ **Hypoallergenic soap**
- ☐ **Fragrance-free moisturizer**
- ☐ **Soft, clean washcloths**
- ☐ **Medical grade adhesive remover, baby oil, or mineral oil**
- ☐ **Bra without wires** (such as a sports bra), if applicable

Keeping your skin healthy may be the best way to prevent skin irritation while using Optune Lua. Proper skin care includes **diligent handwashing, gentle cleansing of the skin, removing hair when appropriate, and extra care when applying and removing adhesive arrays.**

Have questions or need support?



Call MyNovocure at
1-855-281-9301



Email MyNovocure at
support@mynovocure.com

Access resources online



For educational resources about Optune Lua, visit
OptuneLua.com/resources

Please see full Important Safety Information on page 2 and the Optune Lua **Patient Information and Operation Manual (PIOM)** at **OptuneLua.com**.

What is Optune Lua® approved to treat?

Optune Lua is a wearable, portable, FDA-approved device used together with PD-1/PD-L1 inhibitors (immunotherapy) or docetaxel. It is indicated for adult patients with metastatic non-small cell lung cancer (mNSCLC) who have progressed on or after a platinum-based regimen.

Important Safety Information

Who should not use Optune Lua?

Optune Lua for mNSCLC is not for everyone. Talk to your doctor if you have:

- **An electrical implant.** Use of Optune Lua together with electrical implants has not been tested and may cause the implanted device not to work properly
- **A known sensitivity to gels** like the gel used on electrocardiogram (ECG) stickers or transcutaneous electrical nerve stimulation (TENS) electrodes. In this case, skin contact with the gel used with Optune Lua may commonly cause increased redness and itching, and rarely may even lead to severe allergies such as a fall in blood pressure and difficulty breathing

Do not use Optune Lua if you are pregnant or are planning to become pregnant. It is not known if Optune Lua is safe or effective during pregnancy.

What should I know before using Optune Lua?

Optune Lua should only be used after receiving training from qualified personnel, such as your doctor, a nurse, or other medical staff who have completed a training course given by Novocure®, the maker of Optune Lua.

- Do not use any parts that did not come with Optune Lua Treatment Kit sent to you by Novocure or given to you by your doctor
- Do not get the device or transducer arrays wet
- Please be aware that Optune Lua has a cord that plugs into an electrical socket. Be careful of tripping when it's connected
- If you have an underlying serious skin condition where the transducer arrays are placed, discuss with your doctor whether this may prevent or temporarily interfere with Optune Lua treatment

What are the possible side effects of Optune Lua?

The most common side effects of Optune Lua when used together with certain immunotherapy and chemotherapy drugs were dermatitis, pain in the muscles, bones, or joints, fatigue, anemia, alopecia (hair loss), dyspnea, nausea, cough, diarrhea, anorexia, pruritus (itching), leukopenia, pneumonia, respiratory tract infection, localized edema (swelling), rash, pain, constipation, skin ulcers, hypokalemia (low potassium levels), hypoalbuminemia (low albumin levels), hyponatremia (low sodium levels), and dysphagia (difficulty swallowing).

Other potential adverse effects associated with the use of Optune Lua include treatment related skin irritation, allergic reaction to the adhesive or to the gel, overheating of the array leading to pain and/or local skin burns, infections at site where the arrays make contact with the skin, local warmth and tingling sensation beneath the arrays, medical device site reaction, muscle twitching, and skin breakdown/skin ulcer. Talk to your doctor if you have any of these side effects or questions.

Please see Patient Information and Operation Manual (PIOM) for Optune Lua at OptuneLua.com/PIOM.



©2025 Novocure GmbH. All Rights Reserved. Optune Lua, MyNovocure, and Novocure are registered trademarks of Novocure GmbH. US-OPL-00394 v1.0 May 2025

