

Introduction

- Pancreatic cancer remains a leading cause of cancer-related death, with limited treatment options available for the advanced stage.¹
- Tumor Treating Fields (TTFields) are electric fields that have shown *in vitro* and *in vivo* efficacy in pancreatic cancer models, with an optimal frequency of 150 kHz,² and have been shown to induce a state of BRCAness in various cancer types.^{3,4}
- TTFields have improved survival when used together with Gem/NabP in a recent phase 3 clinical trial (PANOVA-3, NCT03377491).⁵
- FOLFIRINOX is a combination regimen of fluorouracil, oxaliplatin, irinotecan and leucovorin, serving as a first-line chemotherapy option in the advanced setting.¹
- BRCA-mutated tumors seem to be more sensitive to FOLFIRINOX.¹
- This *in vitro* study examined the potential concomitant use of TTFields with FOLFIRINOX for pancreatic cancer treatment.

Methods

- TTFields application:** AsPC1 and BxPC3 BRCA-wild type human pancreatic cancer cells were treated with 150 kHz TTFields (1 and 0.7 V/cm RMS, respectively), with or without various FOLFIRINOX concentrations, for 24-72 hrs using the inovitro system.
- Co-application experiments:** TTFields were applied for 72 hrs in the absence or presence of increasing concentrations of FOLFIRINOX, and tested for cell count, colony formation, and apoptosis.
- Cell count:** Treated cells were counted using a flow cytometer, and cell count was calculated relative to control.
- Overall effect:** Treated cells were harvested, re-plated, and grown for an additional 7-14 days. Colonies were stained (0.5% crystal violet solution), counted, and colony formation was calculated relative to control. Overall effect was calculated by multiplying with the corresponding cell count.
- RNA sequencing:** RNA was isolated from TTFields-treated and control samples, and a list of differentially expressed genes (DEGs) by time point and cell type was generated. Gene Set Enrichment Analysis (GSEA) was performed on the gene lists pre-ordered by their z score, using the GSEA software against the Reactome gene set, showing a subset of pathways divided by their Reactome head pathway.
- DNA damage:** Treated cells were fixed (4% paraformaldehyde), permeabilized (0.1% Triton X-100) and stained with anti-γH2AX antibody. Slides were mounted in the presence of DAPI nuclear stain, and images collected on a fluorescent confocal microscope.

References: 1. M.N. Rosen, et al. (2021) *World J Gastroenterol.* 27 (17): 1943-1958. 2. Giladi, M., et al. (2014). *Pancreatology* 14(1): 54-63 3. N.K. Karanam, et al. *Transl. Res.* 2020; 217: 33. 4. H. Mumlat et al. *Lung Cancer.* 2021; 160: 99. 5 H. M. Babiker et al. (2023) *J. Clin. Oncol.* 41(16_suppl): p. TPS4199-TPS4199

Results

FIGURE 1. TTFields augmented the cytotoxic effects of FOLFIRINOX in pancreatic cancer cells




FIGURE 3. TTFields induced DNA damage in pancreatic cancer cells





FIGURE 2. TTFields downregulated cell cycle-related pathways and DNA replication and repair pathways in pancreatic cancer cells



Conclusions

- TTFields application together with FOLFIRINOX in pancreatic cancer cells lacking background BRCA mutations exhibits potential improvement relative to FOLFIRINOX alone.
- The improved efficacy of FOLFIRINOX with TTFields may be rationalized based on the induction of DNA damage and downregulation of key DNA repair pathways.
- Future studies should explore the possibility of reducing FOLFIRINOX treatment dose, potentially alleviating toxicity related to FOLFIRINOX treatment.